


Ackermann<sup>®</sup>  
medical



p|spine

Polyaxial Pedicle Screw System

p|spine  
Table of Content

p spine	Polyaxial Pedicle Screw System	Page
Manufacturer Information	Company Profile	3
	Concentration on the Essentials	3
The p spine Screw System		4
Indications		5
Contraindications		5
Surgical Technique	Patient Positioning and Access	6
	Positioning of the Retractor and Annular Window	6
	Pedicle Opening	7
	Checking of the Pedicle Opening	7
	Determining the Screw Length	8
	Preparing the Pedicle Screw Driver	8
	Valve and Radiator	9
	Screwing in the p spine Screws	9
	Rod Bending	10
	Rod Insertion	10
	Placing the Nut Screw with the Nut Screw Driver	11
	Placing the Nut Screw with the Vertebrae Lever Rocker	11
	Repositioning with the Distraction Forceps	12
	Repositioning with the Compression Forceps	12
	Assembly of the Torque Wrench with T-Handle and Anti-Rotation Forceps	13
	Insertion of the Lumbar Cages	13
Screws	Screws	14
	Cement Screws	14
Instrumentation	Cement Adapter	15
	Crosslinks	15
	Rods Titanium	15
	Rods Nitinol	15
	Instruments	16
	Instrumentation sets	20
Handling & Reprocessing Instructions of reusable Instrumentation		21
Supplementary Equipment		25

# p|spine

Company Profile | Concentration on the Essentials

## Company Profile

We, the Ackermann Medical GmbH & Co. KG, have our head office in Schaffhausen in Switzerland. Since almost one decade we are specialized in the development and production of medical products for spinal surgery. Therefore it is our matter of course to garant our company conforms to the highest medical standards according to FDA, 93/43/EEC (CE), ISO 13485:2003. Since the beginning of 2013 we are building up a direct sales of implants and instruments for spinal surgery in Germany.

## Concentration on the Essentials

We have made it our mission to produce only products offering maximum benefit to patients, the attending physicians, and customers. Therefore we have committed ourselves offering our goods in the best possible quality at consistently low prices.

**This is not only due to our own production and development, but also to our consistent concentration on the essentials**

- no compromises in quality
- optimization of the sales management
- responsible marketing

*Specifications, designs and accessories are subject to change without any notice or obligation on behalf of the manufacturer.*


# p|spine

The p|spine Screw System

Being polyaxial, p|spine pedicle screws allow for an easy positioning and facilitate the fitting. Due to the self-tapping threads pre-drilling becomes obsolete, which later on increases the stability within the vertebra. The specific tulip-shaped design of the self-locking thread prevents the set screw from loosening and provides the possibility of a very low profile. All screws are cannulated and a colour coding ensures that the right size is quickly identified.

The p|spine screws are entirely manufactured from a high grade titanium alloy (TiAl6V4). Rods made of titanium or nitinol support a rigid, as well as a dynamic stabilization of the spinal column.

The p|spine pedicle screws are supplemented by a complete set of instruments, fully cannulated and made in Germany of high quality stainless steel.




# p|spine

Indications | Contraindications

Ackermann provides a full range of sizes with diameters from 5 mm to 7.5 mm and lengths from 25 mm up to 70 mm. Additionally, the p|spine screws are completed by the cement screws, which are available with the diameter 6.5 mm and lengths from 35 mm up to 50 mm.

## Lengths



## Colour-Code and Diameter



## Indications

The use of p|spine polyaxial pedicle screws is intended for:

- degenerative disk disease
- prolapsed intervertebral disc
- pseudarthrosis
- degenerative scoliosis
- revision surgery

## Contraindications

The use of p|spine polyaxial pedicle screws is NOT intended for:

- leukocytosis
- osteoporosis
- pseudarthrosis
- patients with fractures or tumors in the spine area
- patients with spine associated infections
- psychiatric disorder
- pregnancy
- patients with proven material allergy or tendency to react to foreign bodies



## ONE | Patient Positioning and Access

Position the patient in a prone position on an operating table. Use a lumbar support to avoid intraoperative bleeding caused by abdominal compression.

Locate the correct level under x-ray radiation (an x-ray c-arm is recommended) and perform a median incision over the concerned segment. The incision should be done carefully to avoid any case of subcutaneous damage.

### Note

| After dissection, the musculus erector spinae may be separated laterally to obtain the required exposure of the vertebrae and their facet joints.



## TWO | Positioning of the Retractor and Annular Window

After incising and retracting the surrounding tissue, insert the retractor.





For optimal access to the concerned intervertebral disc, perform a laminectomy or laminotomy, and if needed a facetectomy.

Use a nerve root retractor to carefully retract the dura mater and upper nerve roots to the side.

### THREE | Pedicle Opening

First of all the the pedicle walls need to be punctured with a puncture needle or guide wire [70-7781].

Hereby determine the positioning of the screw, which should be inserted in the followin step  
The cannulated awl [70-7796+70-7796B] is to be passed over the fitting puncture needle [70-7796A], by what the cortex is opened..  
The deep rash of the awl prevents the surgeon from penetrating the vertebral body too deep





### FOUR | Checking of the Pedicle Opening

Ensure an optimal screw bearing with the pedicle guide wire [70-7781]. Afterwards exclude a perforation.

### FIVE | Determining the Screw Length

Due to the scale on the depth measuring gauge [70-7794] in combination with 70-7781, the optimal screw length can be determined directly and adequately in mm.




### SIX | Preparing the Pedicle Screw Driver

The pedicle screw driver set is consisting of the following parts (see on the left):

- | trocar for pedicle screw driver [e.g. 70-7780]
- | T-handle [70-7799]
- | and pedicle screw [e.g. 70-7731]


The final assembly looks like the below scheme:






## SEVEN | Valve and Radiator

Although the Ackermann Medical pedicle screws are self-tapping, thread cutters from 5.0 mm to 7.5 mm [70-7757, 70-7758, 70-7783 and 70-7787] are available in the standard product range and also in the operation set. In the case of an extremely hardened bone wall of the cortex, it might be necessary to preturn a thread, to guarantee a firm fitting of the pedicle screws.






## EIGHT | Screwing in the p|spine Screws

All Ackermann Medical pedicle screws are cannulated, so that the screw can be placed over the guide wire into the exact position and can be tightened with sufficient force. The guide wire and the borehole redundantize the axis correction.



### NINE | Rod Bending

If you are not using Ackermann Medical pre-bent connecting rods they can be bent individually with the rod bending forceps [70-7785]



### TEN | Rod Insertion

Insert the (curved) connecting rod [e.g. 70-7771] with the rod holding instrument [70-7795] to restore the lordosis.




### **ELEVEN | Placing the Nut Screw with the Nut Screw Driver**

Prepare the positioning instrument [70-7788+70-7809] for the placement of the nut screw as followst.

Telescope the quick release. Afterwards insert the nut screw. Advance the quick release until it stops. Finally, insert the screw through the guide tube and position it over the rod

#### **Note**

| Do not tighten the screw!



### **Alternatively | Placing the Nut Screw with the Vertebrae Lever Rocker**


Attach the vertebrae lever rocker [70-7708] to the pedicle screw. Then depress the rod with the buttress nut sleeve [see ELEVEN | Placing the Nut Screw with the Nut Screw Driver]. Screw in the nut screw

#### **Note**

| Do not tighten the screw!


**TWELVE** | Repositioning with the  
Distraction Forceps

Use the distraction forceps [70-7797] to distract  
the distance between the two pedicle screws.



**THIRTEEN** | Repositioning with the  
Compression Forceps

Use the compression forceps [70-7798] to  
compress the distance between the two pedicle  
screws.




#### **FOURTEEN** | Assembly of the Torque Wrench with T-Handle and Anti-Rotation Holder

Insert the torque wrench [70-7790] into the anti-rotation holder [70-7791]. Afterwards insert and set the pressure pad on connecting rod. Tightening the nut screw with max. 8 Nm

#### **Note**

| A higher torque can damage the screw! Ackermann Medical screws have a self-locking thread, which automatically connects by 8 Nm with the tulip head.



#### **FIFTEEN** | Insertion of the Lumbar Cages

Now you can start the discectomy as preparation for the insertion of the lumbar cages. For all further steps please see *Surgical Technique t|spine* or *t|spinecurve*.

p|spine  
Screws



Diameter 5.0 mm

Screws	Length [mm]	Colour-Code
70-7720	25	blue
70-7721	30	blue
70-7722	35	blue
70-7723	40	blue
70-7724	45	blue
70-7725	50	blue

Diameter 5.5 mm

Screws	Length [mm]	Colour-Code
70-7730	25	pink
70-7731	30	pink
70-7732	35	pink
70-7733	40	pink
70-7734	45	pink
70-7735	50	pink
70-7737	55	pink
70-7736	60	pink

Diameter 6.5 mm

Schrauben	Length [mm]	Colour-Code
70-7740	30	green
70-7741	35	green



Diameter 6.5 mm

Screws	Length [mm]	Colour-Code
70-7742	40	green
70-7743	45	green
70-7744	50	green
70-7747	55	green
70-7745	60	green

Diameter 7.5 mm

Screws	Length [mm]	Colour-Code
70-7751	40	gold
70-7752	45	gold
70-7753	50	gold
70-7750	55	gold
70-7754	60	gold
70-7755	65	gold
70-7756	70	gold

Diameter 6.5 mm

Cement Screws	Length [mm]	Colour-Code
70-7741CM	35	green
70-7742CM	40	green
70-7743CM	45	green
70-7744CM	50	green

Ackermann offers diverse instruments for save and successfull use of the p|spine screws, which are available separately or in full instrumentation sets.

**Cement Adapter**

70-7749CM      Cement adapter for cement screws

**Crosslinks**

70-7762      Short, 36 mm - 44 mm  
70-7763      Medium, 43 mm - 52 mm  
70-7764      Long, 51 mm - 67 mm

**Diameter 5.5 mm**


**Rods Titanium      Length [mm]**

70-7770      50  
70-7771      60  
70-7772      70  
70-7773      80  
70-7774      90  
70-7775      100  
70-7775-150      150  
70-7775-300      300

**Durchmesser 5.5 mm**

**Rods Nitinol      Length [mm]**

70-7776      60  
70-7777      70  
70-7778      80





### Punction Needle

70-7653 with handle



### Pedicle Finder, 3 mm ball tip, 190 mm WL

70-7687 not malleable



70-7688 malleable



### Awl

70-7695 straight



70-7695SP tapered



70-7694 Forceps for 70-7695/SP



### Nut

70-7719 for crosslinks



70-7718 Nut driver for 70-7719, inner HEX 6 mm for crosslinks



### Stapler

70-7719STAP Stapler for pedicle screw system

### Persuader for Percutaneous Reduction and Rod Replacement

70-7748HE Standard version

70-7748 ► SET PART ECO version



### Tapping Device

70-7787 ► SET PART for 5 mm screws




70-7757 ► SET PART for 5.5 mm screws



70-7758 ► SET PART for 6.5 mm screws



70-7783 for 7.5 mm screws



### Screw Driver for Pedicle Screws

70-7759 ► SET PART Trocar, D 5 mm, for 70-7780

70-7780 ► SET PART Pedicle screw driver, 6 x  
HEX 3.5 mm

70-7781 ► SET PART Guide wire, D 1.3 mm,  
length 400 mm, for 70-7780

### Rod Pushing Device

70-7782 ► SET PART with T-handle

### Rod Bending Plier

70-7785 ► SET PART for Titanium and Nitinol rods

### Screw Driver for Screw Setting

70-7788 ► SET PART

70-7809 ► SET PART Guide sleeve for 70-7788

### Nut Driver

70-7788-100 TX 20

70-7788-200 TX 25

### Torque Wrench

70-7790 ► SET PART with T-handle

### Antirotation Device

70-7791 ► SET PART

### Depth Measuring Gauge

70-7794 ► SET PART

### Rod Holding Device

70-7795 ► SET PART curved

70-7795STR straight

70-7795PCS for percutaneous system



### Nut Driver

70-7788-100 TX 20

70-7788-200 TX 25

### Torque Wrench

70-7790 ► SET PART with T-handle



### Depth Measuring Gauge

70-7794 ► SET PART



### Rod Holding Device

70-7795 ► SET PART curved

70-7795STR straight

70-7795PCS for percutaneous system



### Pedicle Finder

70-7796 ▶ SET PART Trocar, cannulated for 70-7796A/B

70-7796A ▶ SET PART Guide wire for 70-7796

70-7796B ▶ SET PART Holding device for 70-7796



### Distraction Forceps

70-7797 ▶ SET PART



### Compression Forceps

70-7798 ▶ SET PART



### T-Handle

70-7799 ▶ SET PART with ratchet and quick action coupling, cannulated, for 70-7753/57/58/87/88



### Vertebrae Lever Rocker

70-7808 ▶ SET PART



### Instrumentation Sets

►70-7473SET	p spine set with container, insert for instruments and instruments
Assembly	1 x container
	2 x insert for instruments
	3 x tapping device
	1 x rod pushing device
	2 x trocar for pedicle screw driver
	1 x pedicle screw driver
	2 x guide wire
	1 x rod bending forceps
	2 x screw driver for screw setting
	1 x torque wrench
	1 x antirotation device
	1 x depth measuring gauge
	1 x rod holding device
	1 x trocar for guide wire and holding device
	1 x guidance
	1 x holding device
	1 x distraction forceps
	1 x compression forceps
	1 x T-handle
	1 x vertebrae lever rocker

*All components of the kit is autoclavable and reusable. If necessary every component is separately available, too.*



# p|spine

## Handling and Reprocessing Instructions of reusable Instrumentation

### Place of Use (immediately after use)

**IMMEDIATELY** after each use (within no more than 10 min. or before drying of contaminants) the instruments need to be cleaned and impurities removed under running water, using a soft brush or cloth used solely for this purpose. NEVER use a metal brush, steel wool or other cleaning devices containing metal in order to avoid the imminent risk of corrosion. Rinse under cold, running water until all visible impurities and contaminants have been successfully removed.

### Storage and Transport

Place instruments in a container. Keep the inside of the container moist/wet (no contaminants shall dry). Reprocess all instruments soonest possible.

### Preparation for Cleaning

Soak instruments in cold water for at least 5 min. and clean them, using a soft brush or cloth which are being used solely for this purpose. NEVER use a metal brush, steel wool or other cleaning devices containing metal in order to avoid the imminent risk of corrosion. Afterwards, wash down the entire surface of the instrument for 10s. by use of a cleaning gun (min. continuous pressure of 4 bar); articulate moveable parts constantly during the preliminary cleaning. Instruments featuring lumina and/or LuerLock flush channels are to be rinsed for an additional 10s. after visibly clear water has flown from the ports. Place the instruments in an ultrasonic bath for 10min. (35-40kHz for min. 5min. or longer acc. to specifications). Prior to switching on the ultrasonics make sure that all lumina, sheaths, etc. are filled with cleaning fluid!

Note that the preliminary cleaning – even the use of a disinfectant – is only intended as a preparatory step and **DOES NOT** replace the actual disinfection!

### Mechanical Cleaning

Make sure that multiple instruments do not come in contact with each other; especially different materials such as titanium, brass, aluminum, stainless steel, etc. need to be cleaned separately in order to avoid formation of a rust film. Composite instruments (particularly stainless steel combined with ceramics) shall be placed with sufficient distance away from other products to prevent damage arising from the pressure of different thermal expansions.

Instruments have been tested with the following devices:

#### Washer-Disinfector G 7735 CD (Miele)

1. washing cycle: alkaline program (No 105)
2. washing cycle: enzymology program (No 105)

#### Washer-Disinfector G 7836 CD (Miele)

1. two component alkaline/enzymatic program
2. OxiVario

#### Washer-Disinfector Niagara SI PCF (Medisafe) (RECOMMENDED)

1. Cleaning process with pulsed ultrasonic irrigation
2. Cleaning process without pulsed ultrasonic irrigation

The water which is to be used needs to be sterile or nearly sterile (<10 microbes/ml) and low in endotoxins (< 0.25 units/ml). The air which is being used for drying needs to be cleaned by means of micro filters which are regularly checked and maintained. A maintenance schedule has to be documented.

# p|spine

## Handling and Reprocessing Instructions of reusable Instrumentation

### Manual Cleaning

Mechanical cleaning is mandatory with these products as they are classified as class „critical B“ according to the RKI/BfArM-recommendations.

### Disinfection

Take the instruments and place them into the disinfecting bath (Caution: products need to be fully immersed; at least 1cm below the liquid surface).

Multiple instruments shall not come in contact with each other; especially different materials such as titanium, brass, aluminum, stainless steel, etc. need to be disinfected separately in order to avoid formation of a rust film. Composite instruments (particularly stainless steel combined with plastic) need to be disinfected separately in order to prevent damage arising from the pressure of different thermal expansions. Rinse all the lumina of the instrument at least five times using a sterile syringe (min 50ml) and disinfectant.

### After disinfection

Remove products and rinse for at least 5 min. under running water until all disinfectant is removed from the instruments (the water which is to be used needs to be sterile or nearly sterile with <10 microbes/ml and low in endotoxins with < 0.25 units/ml).

Constantly articulate moveable parts.

Rinse all the lumina of the instrument with water at least five times using a sterile syringe (min 50ml).

### Disinfectants that have been successfully tested are

1. Alkaline, Neodisher FA, pH 12.2, Dr. Weigert
2. Enzymatic, deconex 23 Neutrazym, pH 9.7, Borer
3. 2-Component Alkaline/Enzymatic, deconex TWIN PH, pH 10.9, deconex TWINZYME, pH 7, Borer
4. 2-Component Alkaline, Sekumatic FR, pH 12.1; Sekumatic OxiVario. PH 7.8; Neutralizer: Sekumatic FNZ, pH 2.2, Ecolab
5. Enzymatic; M20029 3E-Zyme Scope Plus, pH 6.1, Medisafe
6. Enzymatic + Ultrasound, M20029 3E-Zyme Scope Plus, pH 6.1, Medisafe

### Drying

After cleaning and disinfection place the instruments into suitable containers. Make sure that there is NO residue of the disinfectant.

When drying as part of the cleaning/disinfection cycle is completed make sure that a temperature of 150°C/300°F is not exceeded.

All operations need to take place in a clean, monitored environment!

### Maintenance

Apply a small amount of high-grade surgical lubricant on all joints or other moveable parts which are supposed to move smoothly. Sort out all blunt or damaged instruments.

Clearly damaged instruments (cracks on the insulation, breakage, strongly bleached polymer handles or coatings) are NOT to be reused but repaired or disposed of.

### Testing and Inspection

Jointed instruments are to be tested for ease of movement (avoid too much backlash). The functionality of ratchet mechanisms needs to be checked. All instruments: visually check for damage and wear. Blades should be even and without notches. Long and narrow instruments (especially jointed instruments) should be particularly checked for damages. If instruments are part of a larger set they are to be checked together with all associated components.

# p|spine

## Handling and Reprocessing Instructions of reusable Instrumentation

### Packaging

Individually: a standardized packaging material may be used. The size of each bag needs to match the individual instrument so that there is no tension applied on the sealing.

Sets: sort instruments into designated trays or place on multi-purpose sterilization trays. Blades need to be protected; the weight of each tray may not exceed 8kg (18lbs). For the trays an adequate packaging procedure is to be used.

### Sterilization

All products have been precleaned to an extent which allows for processing and sterilization by use of the equipment stated here. This only applies to a processing method according to these instructions within a system that has been configured and validated in compliance with ISO 17665 and in which all cleaning/disinfecting devices comply with ISO 15883. With the result of the sterilization process greatly depending on the equipment that is being used a sterilization validation acc. to ISO 17665 **MUST** be performed at the place of use prior to the first application. All products **MAY** be used only if the result of this validation is positive.

For the sterilization of medical devices various methods can be applied. Regarding products manufactured by Ackermann, gravity steam sterilization with a fractionated process is recommended.

Temperature	134°C – 137°C (273°F – 279°F)
Pressure	3 bar
Duration	5 min

Please comply with all recommendations issued by the manufacturer of your sterilization device with regard to handling and loading. Instruments that are to be sterilized need to be thoroughly exposed to the steam, including inner surfaces. Before using the instruments they need to be cooled down to room temperature.

Other durations and/or temperatures may also be applied. However, when doing so deviations of parameters should be validated (Caution: contact the manufacturer of your autoclave to confirm temperatures and/or sterilization durations). Temperature inside the autoclave should not exceed 139°C/182°F. This could cause possible damage to handles, insulation or other non-metallic components. Do not sterilize using hot air or Processing and Sterilization Instructions of Medical Devices (acc. to ISO 17664) flash autoclave methods.

In case only pre-vacuum sterilization can be performed, please adhere to the following parameters:

#### for Europe ( except Switzerland and France )

Sterilizer type	pre-vacuum
Preconditioning pulses	3
Preconditioning pressure	30 psia
Minimum temperature	134°C
Cycle time	5 min.
Sample configuration	individually wrapped

#### for Switzerland and France

Sterilizer type	pre-vacuum
Preconditioning pulses	3
Preconditioning pressure	30 psia
Minimum temperature	134°C
Cycle time	18 min.
Sample configuration	individually wrapped

## p|spine

## Handling and Reprocessing Instructions of reusable Instrumentation

## others

Sterilizer type	pre-vacuum
Preconditioning pulses	3
Preconditioning pressure	30 psia
Minimum temperature	132°C / 270 °F
Cycle time	4 min.
Sample configuration	individually wrapped

**Storage**

Store instruments secured against mechanical damage. Use additional wrapping to protect against dust. Do not stack instruments which are packed sterile; especially do not place heavy items on top in order to avoid damage to the sterile packaging of other instruments.

Products need to be stored in a clean and dust-free environment at moderate temperatures of 19° - 25°C (66° - 77°F) and humidity of 40 - 60% (to avoid the risk of embrittlement of the sterile packaging AND of plastic components, especially handles).

**Additional Information**

Do not exceed maximum loading capacity of the sterilizer when processing multiple instruments in one sterilization cycle.

**Warnings****ALL INSTRUMENTS MUST BE CLEANED, DESINFECTED AND STERILIZED PRIOR TO EACH USE.**

All reusable Ackermann products are shipped in non-sterile condition. Before the first use transport packaging, coarse dust/pieces of paper/packaging remains need to be removed and each product processed and sterilized according to these instructions. All products have been pre-cleaned to an extent which allows for processing and sterilization by use of the equipment stated here. This only applies to a processing method according to these instructions within a system that has been configured and validated in compliance with ISO 17665 and in which all cleaning/disinfecting devices comply with ISO 15883.

**THOROUGH CLEANING AND DISINFECTION IS CRUCIAL FOR AN EFFECTIVE STERILIZATION!**

Especially in Germany, the adherence to and knowledge of the RKI/BfArM-recommendations is the prerequisite for these instructions. In accordance with EU Directive 93/42 processing is NOT permitted in Germany without comprehensive awareness of these guidelines!

Processing must only be carried out by personnel explicitly designated by §4 Abs. 3 MPBetreibV after verifying their qualification! Strong cleaning agents may cause fading of markings.

**Limitations on Reprocessing**

Instruments have been validated for 50 cycles, based on an average treatment.

Products that have been marked as single-use (acc. to DIN EN 980 2008-11/figure 5.2) **MUST NOT** be reprocessed! With such products materials are being used that are NOT reprocessable under normal conditions or do not withstand more than one sterilization process and, therefore, may break during surgery if reprocessed! (this applies to practically all single-use products featuring plastic components)



Complementary to our p|spine screws Ackermann offers a wide range of discectomy and neuroscopic instruments.

### Vertibroscope and rotative holding disk

70-7040	Rotative holding disk holds the endoscope and allows 360° rotation with 30° view
70-7180	Vertibroscope for discectomy, 30° rotatable, with light-transmission, lenses made of saphire

### Kerrison rongeur

70-7060	30°, 2 mm jaw, 250 mm WL
70-7065	40°, 2 mm jaw, 250 mm WL
70-7070	90°, 2 mm jaw, 250 mm WL
70-7075	90°, 3.5 mm jaw, 330 mm WL, rotatable
70-7076	40°, 3.5 mm jaw, 330 mm WL, rotatable

### Micro conchotome, 250 mm WL

70-7080	5 mm jaw, 45° upwards arched
70-7085	3 mm jaw, 45° upwards arched
70-7090	5 mm jaw, 45° straight
70-7095	3 mm jaw, 45° straight

### Micro forceps, 2 mm jaw, 200 mm WL

70-7086	45° angled upwards
70-7091	straight
70-7092	straight, spoon with teeth

### Micro scissors, 2 mm jaw, 200 mm WL

70-7097	single action, cutting 90° towards the handle
---------	---

### Curette, 272 mm overall length

70-7100	3 mm jaw, 45° angled
70-7105	5 mm jaw, 45° angled
70-7106	5 mm jaw, 90° angled



Suction tube, 272 mm overall length

70-7110 without spatula tip

70-7111 with spatula tip

Bipolar forceps, 1.2 mm jaw, bayonet, 280 mm overall length

70-7115 straight

70-7120 angled

Bipolar cable

70-7130

Nerve manipulator, 90° angled - knurled handle

70-7135 1 mm tip, 272 overall length

70-7142 3 mm tip, 30 mm spatula length, curved upwards

70-7146 1.5 mm ball tip, 5 mm angled, 200 mm WL

70-7147 1.5 mm ball tip, 9 mm angled, 200 mm WL

Spatula, 90° angled, 200 mm WL - knurled handle

70-7139 spatula size 7 mm x 4 mm, knurled handle

Nerve retractor, knurled handle

70-7140 wide tip, 272 mm overall length


70-7143 3 mm tip, 30 mm spatula length, 200 mm WL

Nerve manipulator, 272 mm overall length

70-7136 1.5 mm tip, 90° angled

70-7137 1 mm tip, 90° angled

70-7141 3 mm jaw, spatula form



Ball tip dissector, 272 mm overall length

70-7145                      90°, 1.5 mm tip

Bayonet knife, 272 mm overall length

70-7150                      2 mm tip

70-7155                      1 mm tip

Neurosurgical retraction arm

70-7190                      to hold the rotative disk and the  
vertebroscope



Ackermann Medical GmbH & Co. KG  
Schlossstrasse 67  
8200 Schaffhausen  
Switzerland

Tel +41 (0) 525 330 915  
Fax +41 (0) 525 330 948  
info@ackermannmedical.ch  
www.ackermannmedical.com