

CATALOGUE 2022

ACCESSORIES



ACCESSORIES

MODULAR AUTO LOCK MECHANISM

What Is Auto Lock Mechanism

Auto lock mechanism is a heat-sensitive automatic locking system developed in order to ease sterilization control by R&D department. Modular auto lock system designed to be usable with different container systems.

This system has been developed in accordance with the requirements of sterilization to be used in a steam autoclave for our container systems. Auto lock is a system designed in reusable property instead of disposable plastic products that are used as security locks in containers. The system has been tested in our own facilities 5000 cycles and has successfully passed al tests. Minimum expected life of the system is 5000 sterilizations. During this time system is usable without the need for any maintenance or service.



How does Auto Lock Mechanism Work?

The average temperature in the autoclave during sterilization is at the level of 135 $^\circ$ C degrees.


As the temperature in the autoclave begins to rise above 100 degrees, automatic locking mechanism located on the box is activated. When the mechanism is activated, the red movable locking shaft in it is going forward (Picture 1). When the movable locking shaft is going forward, it blocks the locking mechanism in lid and so it prevents uncontrolled opening of the lid.

The process of the system becoming inactive is as following: After sterilization process has been completed and container has been removed from the autoclave, the container and Auto Lock system are allowed to cool between 5 and 10 minutes. Then the red movable locking shaft is pushed back with the help of fingers (Picture 3). Thus, lid closure of the container ceases to be blocked and comes to opened position.

When the container enters into each sterilization, these processes are repeated in the same manner.

When the user sees that Auto Lock mechanism is activated (Picture 1), he/she can understand that the container has entered to the sterilization process. If Auto lock mechanism is inactive (Picture 2), it means that this container has been previously opened.

Picture 1 (Closed Position)



Picture 2

WARNING: It's possible to move the movable shaft with only the effect of temperature. It's not possible under normal conditions to activate (to be brought into locked position) lock mechanism.

What Are The Advantages Of Modular Auto Mechanism?

Because of system being reusable, it is more economical when you compare with plastic security seal usage.


System is modular, so it can be usable with different brand products. In this reason, so many users can make use of system advantages.

After the usage of disposable products, plastic material harms the nature. Because of reusable property of auto lock mechanism and because of recycling raw material inside, system does not harm the nature.

Disposable products have many disadvantages such as purchasing, storage and inventory tracking. Auto Lock mechanism saves the user from this workload and cost.

Disposable products must be attached to the (2) container(s) prior to each sterilization and this leads to loss of time for the employees.

Auto Lock Mechanism removes the difficulties faced during the use of security seals. In breakdown condition, because of product is not fixed on the container, system can be easily serviceable or changeable so container does not needs replacing.




Picture 4



All our container systems with Auto Lock mechanism are designed to be used also with plastic security seal when necessary (Picture 4).

MODULAR AUTO LOCK MECHANISM



Picture 1



Locking clamp should be removed from the groove before assembly.

Picture 3



When you sure that product is snug, locking clamp should be placed through the arrow direction.



Assembly And Disassembly Of Modular **Auto Lock Mechanism**

Picture 2



Then it should be placed to the container label part as shown picture.

Picture 4



Please be sure that locking stamp is forwarded to the end of the clamp groove and it is locked.

Lock mechanism should be in open position for the disassembly process (Picture 2).

<u>It is not possible</u> to make disassembly on closed position (Picture 1).

When you sure that lock mechanism is on open position, locking clamp can be removable by pulling backward through the reverse direction of the locking direction.

STERILE CONTAINER SYSTEMS

WINDOW FOR **EMULATING INDICATOR**

We continue to produce solution for the sterilization management requirements with our 20 years of experience in medical industry. For this purpose, we have developed emulating indicator window for the sterile container systems.

Aim Of The Window For Emulating Indicator

Emulating indicators are used for the load sterilization monitorization while surgical instruments are loaded to containers in central sterilization services department (CSSD).

During standard container usage it is not possible to see indicator results without opening the lid. When set is decided to use in surgery if result is negative and also if there is no alternative set, surgery is delaying and problems occurs in implementations.

Result of the emulating indicator can observed after cycle easily from the window on the container case. In negative conditions, sterilization processes are repeated before send set to usage. By this way waste of time prohibited and it can help to avoid troubles before start.

Window For Emulating Indicator

- Transparent PPSU polymer material is used on window part of the system which is suitable to autoclave. PPSU polymer material has resistance to scoring and impact also it is a special material which has heat resistance up to 330°
- Silicon joint is used between window and box to avoid sterilization risk inside the container and security is provided %100 by 4 pieces rivet connection.
- All other parts used in the system is manufactured from anodized aluminum and stainless steel.
- Indicator is designed to place the box from two directions. Also there is a spring mechanism to immobilize indicator because of the pressure inside the indicator placed sections. This spring mechanism and the immobilization property from the two side of indicator are provided to use different branded indicators to use in this system together.











Usage Of Window For Emulating Indicator

Emulating indicator before placing materials to sterilize *inside the box for the sterilization purpose is placed inside* the box as below.

When you placed indicator band inside the box, you should be sure that placed references and indicator should appears. You should specify the suitable position by reciprocating motions or rotating reverse the indicator band.

After you placed indicator band fairly, you can start sterilization process by placing materials to sterilize after closing the lid.

ACCESSORIES



STR NO STR DATE DUE DATE STR DATE DUE DATE NAME -® • 2 STEAM BLACK NAME SG 1000 ISO EN 11140-1 **Code** CN-2673 **Code** CN-2672 100 pcs. With indicator. 100 pcs. With indicator. For full and Half size containers For small size and special containers Identification Labels (Aluminium) Color Code Thoracotomy Red Blue Basic Set Thoracotomy Green Yellow CN-2703 Grey Thoracotomy Basic Set Black Purple CN-2706 Basic Set Bronze CN-2707

Label (Papers)

Coloured identification labels Engraved to your specifications Maximum 16 letters

Identification (Aluminium)

Grey coloured recommended for use in Plasma.Color fading By effect of Hydro-gen Peroxide Gas Plasma on other colors.






CN-2700 CN-2701 CN-2702 CN-2704 CN-2705



LABELS

Color	,	Code
	Red	CN-2570
	Blue	CN-2571
	Green	CN-2572
	Yellow	CN-2573
	Grey	CN-2574
	Black	CN-2575
	Purple	CN-2576
	Bronze	CN-2577

Coloured identification labels Engraved to your specification








plasma sterilization



For Dental Model	Paper Filter	Paper Filter	Teflon Filter	FILTER RETAINER
	rangen Barangenten Karangenten Karang Barang	•		Code CN-2720 Diameter 186 mm Filter Retainer, for 1/1, 3/4, 1/2, Extra large, Extra long, Large dental
	Code CN-3612 Dimension 165 x 235 mm Single use 100 pcs with indicator.	Code CN-3615 Dimension 165 x 235 mm Single use 100 pcs with indicator. No marking.	CodeCN-2743Dimension 150 x 233 mmNumber of sterilization cycle > 2000	Code CN-2722 Diameter 186 mm Filter Retainer, for all A1 models








Code CN-21716

Dimension 243 x 120 mm

Filter Retainer, for Mini model



Filter Retainer, for 1/2 Dental model




Code CN-2730

Dimension 200 x 200 x 64 mm

Filter Dispenser



Filter Dispenser



Code CN-21710 Dimension 169 x 151 x 53 mm

Filter Dispenser






Code	Dimension
CN 2510	460 405
CN-2519	160 x 135 mm
CN-2520	275 x 125 mm
CN-2521	250 x 235 mm
CN-2522	340 x 230 mm (2 pcs)
CN-2523	300 x 300 mm
CN-2524	380 x 230 mm

Code	Dimension
CN-2525	440 x 230 mm
CN-2526	520 x 230 mm
CN-2527	230 x 220 mm
CN-2528	475 x 230 mm
CN-2529	260 x 180 mm
CN-2532	280 x 200 mm





Code	Dimension
CN-2530	250 x 240 mm
CN-2531	485 x 250 mm
CN-2533	534 x 247 mm
CN-2534	395 x 245 mm
CN-2535	527 x 335 mm
CN-2536	658 x 247 mm

Silicone Mesh Mat, Enables surgical instruments to be stored without slip and dry well.





SILICONE GASKET



Code	Model	Size
CN-3948	Large Dental Silicone Gasket	118 cm
CN-3949	1/2 Dental Silicone Gasket	72 cm
CN-3950	Mini Silicone Gasket	77 cm
CN-3951	Dental Silicone Gasket	89 cm
CN-3952	Flat Silicone Gasket	101 cm
CN-3953	1/2 Silicone Gasket	101 cm
CN-3954	3/4 Silicone Gasket	136 cm
CN-3955	1/1 Silicone Gasket	158 cm
CN-3956	Extra Large Silicone Gasket	180 cm
CN-3957	Extra Long Silicone Gasket	185 cm
CN-3958	Endo Silicone Gasket	130 cm





Code	Outside Dimension
CN-31090	502 x 255 x 130 mm

Rack, for endoscopy

• Wire basket 1/1 50 mm

• Container 1/1 200 mm

• Silicone mat 520 x 230 mm

• Holds 10 endoscopic instruments

• Fits into a 1/1 size 20 cm deep container

Code	Outside Dimension	
CN-31093	555 x 255 x 130 mm	

Rack, for endoscopy

- Holds 10 endoscopic instruments
- Fits into a Extra Long size 18 cm deep container
- Wire basket Extra Long 50 mm
- Container Extra Long 180 mm
- Silicone mat 340x230 (2pcs)

Code	Outside Dimension	
CN-31092	515 x 250 x 165 mm	
Rack, for artrosco	ру	
• Holds 20 artros	copy punches	
• Fits into a 1/1 size 20 cm deep container		
• Wire basket 1/1 50 mm		

- Container 1/1 200 mm
- Silicone mat 440 x 230 mm













Code	Outside Dimension
CN-31094	220 x 250 x 160 mm

Rack, for kerrison • Holds 10 kerrison punches • Fits into a 1/2, 1/1 20 cm and Extra Long size 18 cm deep container • Wire basket 1/2 50 mm • Wire basket 1/1 50 mm

- Wire basket Extra long 50 mm
- Container 1/2 200 mm
- Container 1/1 200 mm
- Container Extra long 180 mm

RACKS FOR ENDOSCOPY














Code

CN-1680-02

CN-1680-03 (1 Piece)

(1 Piece)

٩

Kep nut

A state

Screw

Actual Size

			-
	Code	Lenght	
	CN-12106	238 mm	2
With rail, screw and nut	CN-12107	150 mm	4
	CN-12108	119 mm	4
	CN-12147	87 mm	2



	Code	Lenght
	CN-11601	250 x 9 mm
_	CN-11602	158 x 9 mm
	CN-11603	123 x 9 mm
_	CN-11604	97 x 9 mm
	CN-11605	73 x 9 mm
	CN-11606	23,5 x 9 mm

Can be manufactured in any size demanded



With rail,	screw and nut	

With rail, screw and nut

CN-12111 112 mm

CN-12148 77 mm

Code

CN-12109

CN-12110

Code	Lenght	Height
CN-12100	250 mm	45 mm
CN-12101	152 mm	45 mm
CN-12102	125 mm	45 mm
CN-12145	97 mm	45 mm

Lenght

238 mm

158 mm

Height

45 mm


45 mm

45 mm

45 mm







	Code	Lenght
	CN-12103	240 mm
With rail, screw and nut	CN-12104	150 mm
	CN-12105	120 mm
	CN-12146	90 mm

Instruments holder, from Ø 3 to 8 mm With rail, screw and nut





Code	Lenght
CN-12113	73 mm


Instruments holder, Ø 5 mm With rail, screw and nut

|--|


SILICONE FIXATION SYSTEM RAILS











SILICONE FIXATION SYSTEM RAILS



		Actual Size		
	CodeLenghtHeightCN-1213072 mm34 mmEndoscope holder, doubleOne side Ø 12 mm openOne side from Ø 3 to Ø 15 mmOne side from Ø 3 to Ø 15 mm	Ø 12 mm		CodeLenghtHCN-1213538 mm33Endoscope holder, Ø 4 mm clWith rail, screw and nut
	Code Lought Hitch	Actual Size		
*	CodeLenghtHeightCN-1214172 mm34 mm		•	CodeLenghtHCN-1213638 mm33
	Endoscope holder, double One side Ø 12 mm open One side from Ø 3 to Ø 8 mm closed With rail, screw and nut	Ø 12 mm		Endoscope holder, Ø 6 mm cl With rail, screw and nut
		Actual Size		
	Code Lenght Height			CodeLenghtHCN-1213738 mm33
	CN-12131 72 mm 34 mm Instruments holder from Ø 3 to Ø 15 mm With rail, screw and nut			Endoscope holder, Ø 8 mm cl With rail, screw and nut
		Actual Size		
* * * *	CodeLenghtHeightCN-1213272 mm24 mm	23232323		CodeLenghtHCN-1213838 mm33
	Instruments holder, from Ø 3 to Ø 8 mm With rail, screw and nut			Endoscope holder, Ø 12 mm With rail, screw and nut
		Actual Size		
	CodeLenghtHeightCN-1213338 mm24 mm			Code Lenght H
	Endoscope holder, double one side	Ens Ens		CN-12139 38 mm 33
	from Ø 3 to Ø 8 mm With rail, screw and nut			Endoscope holder, Ø 12 mm With rail, screw and nut
	Code Lenght Height	Actual Size		
	CN-12134 28 mm 34 mm			Code Lenght H CN-12140 38 mm 33
	Endoscope holder, single one side Ø			
Y	12 mm open one side from Ø 3 to Ø 15 mm With rail, screw and nut			Endoscope holder, Ø 15 mm With rail, screw and nut

SILICONE FIXATION SYSTEM RAILS







FIXATION CLAMP

All rights reserved.No part of this catalog may be reproduced in any from without the written permission from ASELTECH

CENTRAL FACTORY

Köprübaşı Mahallesi Selanik Caddesi No:255 Tekkeköy/Samsun/Turkey

T: +90 362 266 59 29 F: +90 362 266 59 39

> info@aseltech.com.tr



www.**aseltech**.com.tr