

# Classic Portable Autoclave



 **Prestige Medical**®  
Integrated decontamination solutions



Prestige Medical Classic – with more than 1 million sold in 80 countries worldwide, you can rely on our reputation for quality and performance.

## Technical Specification

Classic Portable Autoclave							
Model Number	210001	210002	210004	210047	210048	210134	
Max width - overall	340mm	340mm	340mm	340mm	340mm	340mm	340mm
Max height - overall	335mm	420mm	420mm	335mm	420mm	420mm	335mm
Max net weight	4.0kg	4.8kg	4.8kg	4.0kg	4.8kg	4.8kg	5.4kg
Chamber diameter	210mm	210mm	210mm	210mm	210mm	210mm	210mm
Chamber height	236mm	328mm	328mm	236mm	328mm	328mm	236mm
Chamber capacity	9 litres	12 litres	12 litres	9 litres	12 litres	12 litres	9 litres
Max instrument length	228mm	290mm	290mm	228mm	290mm	290mm	228mm
Max load weight	3kg	4kg	4kg	3kg	4kg	4kg	3kg
Sterilizing temperature	126°C	126°C	126°C	121°C	121°C	121°C	134°C
Min sterilizing time	11 minutes	11 minutes	11 minutes	28 minutes	28 minutes	28 minutes	3 minutes
Cycle time - normal	22 minutes	22 minutes	22 minutes	70 minutes	70 minutes	70 minutes	20 minutes
Operating pressure	1.4 bar	1.4 bar	1.4 bar	1.05 bar	1.05 bar	1.05 bar	2.2 bar
Voltage	230v	230v	230v	230v	230v	230v	230v
Wattage	1250w	1250w	1250w	1250w	1250w	1250w	1250w
Frequency	50-60Hz	50-60Hz	50-60Hz	50-60Hz	50-60Hz	50-60Hz	50-60Hz
Printer							
Temperature & Pressure Gauges			✓				✓
Furniture	219708	219720	219720	219708	219720	219720	219708

All Classic units are supplied complete with electrical cord set and instruction manual as appropriate.



Cassette x 3 &  
Carrier - ref: 219706

Cassette  
ref: 219291

Cassette  
Carrier  
ref: 219740



Cord Sets  
UK ref: 219258  
Euro ref: 219297  
UL ref: 219299

**A wide choice of accessories** - allows you to specify the model to suit your exact needs.

Standard  
Basket  
ref: 219708



**Reliable and trusted worldwide**

Designed to give you the most cost effective sterilization solution with a wide choice of models and accessories that are guaranteed to meet your needs.

Steam sterilization, or autoclaving, is the preferred method of eliminating the risk of cross infection through contaminated instruments. It is the only method recommended throughout the world by healthcare professionals including the World Health Organisation (WHO) and the Centre for Disease Control (CDC).

Sterilizing Indicator Strips



126°C ref: 219277  
121°C/134°C ref: 259277

Extended  
Basket  
ref: 219720



Lifter Fork  
ref: 219294



Green Gaskets x 1  
ref: 219500



## Prestige Medical - integrated decontamination solutions

Prestige Medical Ltd is an ISO 9001-2000/EN13485 quality accredited company. The Classic is designed to the latest international standards for safety and performance and is fully CE marked to the Medical Devices Directive 93/42/EEC.



---

### Contact us for advice on infection control



Prestige Medical Limited  
East House, Duttons Way,  
Shadsworth Business Park,  
Blackburn BB1 2QR

Tel : +44 (0) 1254 682 622  
Fax : +44 (0) 1254 682 606

[www.prestigemedical.co.uk](http://www.prestigemedical.co.uk)  
[sales@prestigemedical.co.uk](mailto:sales@prestigemedical.co.uk)

 @PrestigeMedica1

Registered in England  
Reg No. 2826793  
Vat No. 633707835

Registered Office:  
Barlow Road,  
Aldermans Green Industrial Estate,  
Coventry CV2 2LD

Prestige Medical is a registered trademark

Due to a continuous process of development, Prestige Medical reserves the right to amend the specification of products and services without notice