



Prestige Medical



Ultrasonic Cleaner

The Prestige Medical Ultrasonic Cleaner is ideal for the cleaning of a broad range of instruments prior to washing and sterilisation. It is also suitable for using within enzymatic situations. For more information please contact the head office.

Technical Specification

NAME	Ultrasonic Cleaner	POWER SUPPLY	<ul style="list-style-type: none">• AC100-120V60Hz• AC200-240V 50Hz
MODEL	4820	POWER	70W
OVERALL SIZE	290x223x185mm	FREQUENCY	40KHz
TANK SIZE	250x150x80mm	TIME SETTING	0-8Min
VOLUME	2500ml	G.W	2538±10g
Max Temp	42.5°C		

- Fast Easy Operation
- Intuitive Display
- Lightweight
- 2500ml Capacity
- Warming Function to 42.5°C

Get in touch:

www.prestigemedical.co.uk

sales@prestigemedical.co.uk

Tel: +44 (0) 1254 682 622

East House, Duttons Way,

Blackburn, England BB1 2QR





Prestige Medical

- Fast & Practical Operation
- Compact Design
- Automatic Heating
- Seal Length 250mm
- Seal Width 12mm
- Wall Mountable



Pouch Sealer

The Prestige Medical Pouch Sealer is designed to ensure air and water tight packaging of loads ready for sterilisation. For more information on our Pouch Sealer please contact the head office.

Technical Specification

Packet tray		Optional	Dimensions		
Electrical supply voltage and frequency	V-Hz	Volt 230 – Hz 50/60	Width	mm	470
Power	W	100	Height	mm	300
Length of seal	mm	250	Depth	mm	270
Width of seal	mm	12	Packaging dimensions:		
Heating		Pre-set	Width	mm	480
Start-up		Manual	Height	mm	310
Weight	kg	6	Depth	mm	280

Get in touch:

www.prestigemedical.co.uk

sales@prestigemedical.co.uk

Tel: +44 (0) 1254 682 622

East House, Duttons Way,

Blackburn, England BB1 2QR

