

Finest Surgical Technology Since 1954

Suture Passer

The Ackermann Suture Passer is a fascia closing device with stop function that offers a procedural cost effective solution and allows for safe and easy peritoneal closure. This device provides a full thickness closure by minimizing incisional hernia. Suitable for standard and bariatric surgeries, this fascial closure technique provides physicians with a quick and reliable solution in laparoscopic surgeries.



Suture Passer



SUTURE Features &



KEY FEATURES

- ✓ Medical grade – High quality stainless steel
- ✓ 100% made in Germany
- ✓ Simple but yet smart and lightweight engineering for optimal comfort and ease of use
- ✓ Working length of 180mm
- ✓ Small diameter of 2,8mm – for a minimally invasive approach
- ✓ Needle cone shaped tip to facilitate insertion
- ✓ Ergonomic spring handle suitable for left and right-handed physicians
- ✓ Luer lock connector - Flusing port of 45° for easy reprocessing
- ✓ Single action serrated inner jaws for optimal ligature grip
- ✓ Large opening of movable tip to facilitate gripping of the suture





PASSER & Benefits



KEY BENEFITS

- ✓ Cost effective solution
- ✓ Can be used with various laparoscopic suture USP (Unites States Pharmakopeia)
- ✓ Quick and easy suturing of trocar openings and thick fascial layers
- ✓ Full wound closure without removal of trocars
- ✓ Peritoneum is maintained due to small diameter
- ✓ The stop function allows for a controlled insertion and prevents injuries to tissues or organs
- ✓ Reduces risk of hernia and complications
- ✓ Suitable for standard and bariatric procedures
- ✓ Minimally invasive solution for an even quicker healing process



COMPLEMENTARY INSTRUMENT (OPTIONAL)

ITEM # 13-1591ET

A 5 mm straight needle holder with ratchet which is used to temporarily grasp the suture inside the peritoneal cavity after the first insertion by the Fascia Closing Device.