



Urology

Cystoscopy, Resectoscopy & other Specialities

Our Vision

"True to our vision **INNOVATING TRADITION**, it has always been our goal to not only provide our customers and healthcare professionals with an industry-leading value proposition through continuous innovation in our products and processes.


Ackermann

We believe that it is equally important to not lose sight of our roots and to preserve the craftsmanship that has been perfected over generations as an established family business. This enables us to maintain a productive and sustainable home for our staff and to improve clinical outcomes every single day in over 100 countries around the world."

Peter Grassl, CEO

+ 20k SKUs

As **Full-Range-System Manufacturer** we are able to provide you with a coherent product portfolio - from endoscopy to ligation, from implantology to general surgery.



Capabilities

+ Manufacturing


Surgical Instruments, Units, Implants, Optics

+ Regulatory

ISO 13485:2016, US FDA, CFDA, ANVISA,
COFEPRIS


+ Service

+ R&D



Global Player

Sales subsidiaries located in
Greece, India, USA & Mexico



Serving over
100
Countries

More than
80
Employees



3
Manufacturing
Sites

400
Strategic
Partners



Table of Contents

Product Range at a Glance



Urology Sets

Cystoscopy Sets

17 / 23 Fr. Cystoscopy Set	p. 6
Mini Compact Cystoscopy Set	p. 7

Resectoscopy Sets

24 / 26 Fr. Bipolar Resectoscopy Set	p. 8
24 / 26 Fr. Universal Resectoscopy Set	p. 9
17.5 / 18.5 Fr. Universal Resectoscopy Set	p. 10
27 / 28.5 Fr. Monopolar Resectos Set	p. 11
24 / 26 Fr. Monopolar Resecto Set	p. 12
19 / 22 Fr. Monopolar Resecto Set	p. 13

Laser Surgery Instrumentation Sets

23 Fr. Universal Laser Cysto Set	p. 14
24 / 26 Fr. Laser Resecto Set	p. 15
24 / 26 Fr. Laser Resecto Set, demountable	p. 16
24 / 26 Fr. Laser Resecto Set, 0.8 mm	p. 17
24 / 26 Fr. Laser Resecto Set, 1.2 mm	p. 18
19 / 22 Fr. Laser Resecto Set, 0.8 mm	p. 19
19 / 22 Fr. Laser Resecto Set, 1.2 mm	p. 20
17.5 / 18.5 Fr. Laser Resecto Set, 0.8 mm	p. 21
17.5 / 18.5 Fr. Laser Resecto Set, 1.2 mm	p. 22

Lithotripsy Sets

Mechanical Lithotripsy Set	p. 23
----------------------------	-------

Urethrotomy Sets

21 Fr. Visual Urethrotomy Set	p. 24
15 Fr. Visual Urethrotomy Set	p. 25

Nephroscopy Sets

24 / 26 Fr. Short Nephroscopy Set	p. 26
24 / 26 Fr. Standard Nephroscopy Set	p. 27
24 / 26 Fr. Compact Nephroscopy Set	p. 28
24 / 26 Fr. Mini Nephroscopy Set	p. 29

Uretero-Renoscopy Sets

FullHD Uretero-Renoscopy Set	p. 30
Thin HD Uretero-Renoscopy Set	p. 31
Ultra-Thin HD Uretero-Renoscopy Set	p. 32

Cystoscopy

Adults & Adolescents

FullHD Cystoscopes	p. 36
Cysto-Uretero Sheaths & Obturators	p. 38
Optical Instruments	p. 41
Telescope Bridges	p. 43
Albarran Deflectors	p. 44
Flexible & Semi-Rigid Instruments	p. 46
Flexible Monopolar Electrodes	p. 48
Flexible Bipolar Electrodes	p. 49
Dormia Stone Baskets	p. 49

Pediatrics

Compact Mini Cystoscopes	p. 50
Cysto-Uretero Sheaths & Obturators	p. 52
Flexible & Semi-Rigid Instruments	p. 53
Monopolar Electrodes	p. 54
Flexible Cystoscope	p. 57

Resectoscopy

Bipolar Resectoscopy

Bipolar Working Elements	p. 61
Bipolar Electrodes	p. 61
Resectoscopy Sheaths & Obturators	p. 63

Universal Resectoscopy

Universal Working Elements	p. 69
Universal Electrodes	p. 70
Resectoscopy Sheaths & Obturators	p. 72

Monopolar Resectoscopy

Monopolar Working Elements	p. 77
Monopolar Electrodes	p. 77
Resectoscopy Sheaths & Obturators	p. 80

Lithotripsy

Lithotripsy Sheaths & Obturators	p. 90
Mauermayer Stone Punch	p. 91
Stone Lithotripter "The Beast"	p. 91
Optical Stone Forceps & Adapter	p. 91

Urethrotomy

Adults & Adolescents

Working Elements	p. 96
Urethrotomy Sheaths & Obturators	p. 97
Cold Knives	p. 99
Otis Urethrotome	p. 100

Pediatrics

Working Elements	p. 101
Urethrotomy Sheaths & Obturators	p. 101
Cold Knives	p. 102

Nephroscopy

Short & Standard Nephroscopy

Short & Standard Nephroscopes	p. 106
Nephroscopy Sheaths & Obturators	p. 108
Rigid Instruments	p. 109

Compact Nephroscopy

Compact Nephroscope	p. 112
Compact Nephro. Sheath & Obturator	p. 113
Rigid Instruments	p. 114

Mini Nephroscopy

Mini Nephroscope	p. 116
Mini Nephro. Sheaths & Obturators	p. 117
Semi-Rigid Instruments	p. 118
Dilators & Guide Wires	p. 120

Uretero-Renoscapy

Uretero-Renoscopes	p. 124
Flexible & Semi-Rigid Instruments	p. 126
Flexible Uretero-Renoscope	p. 131

Accessories for Urology

Speciality Instruments


Bougies by Dittel	p. 134
Bladder Syringe	p. 135
Strauss Penis Clamp	p. 136
Ellik Evacuator	p. 136

Connection Cables

Bipolar Cables	p. 137
Universal Cable	p. 137
Monopolar Cables	p. 138

17 / 23 Fr. Cystoscopy Set

- | | |
|---|---------------|
| 1x FullHD Cystoscope, Ø 4 mm, 30°, 300 mm | # 32-4020HD |
| 1x Cysto sheath incl. obturator, 2 lateral stopcocks, Ø 17 Fr. | # 32-4070 |
| 1x Cysto sheath incl. obturator, 2 lateral stopcocks, Ø 23 Fr. | # 32-4085 |
| 1x Telescope bridge with 2 instrument channels | # 32-4105 |
| 1x Albarran deflector with 2 instrument channels, lockable | # 32-4115 |
| 1x Flexible ball electrode, monopolar, Ø 5 Fr., 470 mm | # 32-4202 |
| 1x Flexible ball electrode, monopolar, Ø 7 Fr., 470 mm | # 32-4150 |
| 1x Biopsy forceps, Ø 7 Fr., flexible, 400 mm | # 32-4120 |
| 1x Serrated grasping forceps, Ø 7 Fr., flexible, 400 mm | # 32-4130 |
| 1x Scissors, Ø 7 Fr., flexible, 400 mm | # 32-4136 |
| 1x Dormia stone basket, Ø 5 Fr., 4 wires, 600 mm | # 32-4204 |
| 1x Handle for dormia stone basket, Ø 5 Fr | # 32-4211 |
| 1x Optical biopsy spoon, 4.3 mm, double-action, spoon-type | # 32-4055 |
| 1x Bladder syringe, 150 ml | # 32-4395-150 |
| 1x Ellik evacuator | # 32-4106 |
| 1x Monopolar cable, for use with monopolar Urology instruments, 3 m | # 11-1260UR |



Mini Compact Cystoscopy Set

For Pediatrics


- | | |
|--|-----------|
| 1x Mini Compact Cystoscope, for instruments up to Ø 4 mm | # 32-4012 |
| 1x Serrated grasping forceps, Ø 3 Fr., 280 mm | # 32-4160 |
| 1x Biopsy forceps, Ø 3 Fr., 280 mm | # 32-4162 |



24 / 26 Fr. Bipolar Resectoscopy Set

For Adults


- 1x FullHD Cystoscope, Ø 4 mm, 30°, 300 mm # 32-4020HD
- 1x Resectoscopy sheath incl. obturator, continuous-flow, rotatable, Ø 24 / 26 Fr. # 32-4205RR
- 1x Visual obturator with Ø 8 Fr. working channel for Ø 24 / 26 Fr. resectoscopy sheath # 32-4228-6
- 1x Bipolar working element, closed handle, passive # 32-4220BP
- 1x Loop electrode 30°, Ø 24 Fr., bipolar # 32-4240BP
- 1x Ball electrode Ø 3 mm, Ø 24 Fr., bipolar # 32-4250BP
- 1x Knife electrode , Ø 24 Fr., bipolar # 32-4265BP
- 1x Protection- and autoclaving tube for electrodes and cold knives # 32-4230
- 1x Bipolar cable for use with bipolar resectoscope, 3 m # 13-1641UR



24 / 26 Fr. Universal Resectoscopy Set

For Adults


- | | |
|--|----------------|
| 1x FullHD Cystoscope, Ø 4 mm, 30°, 300 mm | # 32-4020HD |
| 1x Resectoscopy sheath incl. obturator, continuous-flow, rotatable, Ø 24 / 26 Fr. | # 32-4205RR |
| 1x Visual obturator with Ø 8 Fr. working channel for Ø 24 / 26 Fr. resectoscopy sheath | # 32-4228-6 |
| 1x Universal working element, closed handle, passive | # 32-4220U |
| 1x Loop electrode 30°, Ø 24 Fr., universal | # 32-4240U |
| 1x Ball electrode Ø 3 mm, Ø 24 Fr., universal | # 32-4250U |
| 1x Knife electrode , Ø 24 Fr., universal | # 32-4265U |
| 1x Protection- and autoclaving tube for electrodes and cold knives | # 32-4230 |
| 1x Neutral rubber electrode for adults, round plug, 5 m | # 16-2000-1212 |
| 1x Universal cable for use with universal resectoscope, 3 m | # 13-1641U |



17.5 / 18.5 Fr. Universal Resectoscopy Set

For Adolescents


- 1 x FullHD Cystoscope, Ø 2.9 mm, 30°, 300 mm # 32-4022HD
- 1 x Resectoscopy sheath incl. obturator, continuous-flow, rotatable, Ø 17.5 / 18.5 Fr., QL # 32-4282RR
- 1 x Visual obturator with Ø 5 Fr. working channel for Ø 17.5 / 18.5 Fr. resectoscopy sheath, QL # 32-4342-17
- 1 x Universal working element, closed handle, passive, QL # 32-4275U
- 1 x Loop electrode 30°, Ø 17.5 Fr., universal # 32-4287U
- 1 x Ball electrode Ø 3 mm, Ø 17.5 Fr., universal # 32-4288U
- 1 x Knife electrode, Ø 17.5 Fr., universal # 32-4289U
- 1 x Protection- and autoclaving tube for electrodes and cold knives # 32-4230
- 1 x Neutral rubber electrode for adults, round plug, 5 m # 16-2000-1212
- 1 x Universal cable for use with universal resectoscope, 3 m # 13-1641U



27 / 28.5 Fr. Monopolar Resectoscopy Set

For Adults


- | | |
|--|---------------|
| 1x FullHD Cystoscope, Ø 4 mm, 30°, 300 mm | # 32-4020HD |
| 1x Resectoscopy sheath incl. obturator, continuous-flow, rotatable, Ø 27 / 28.5 Fr. | # 32-4221RR |
| 1x Visual obturator with Ø 9 Fr. working channel for Ø 27 / 28.5 Fr. resectoscopy sheath | # 32-4240-9 |
| 1x Monopolar working element, closed handle | # 32-4220 |
| 1x Loop electrode 90°, Ø 27 Fr., monopolar | # 32-4295 |
| 1x Loop electrode 30°, Ø 27 Fr., monopolar | # 32-4300 |
| 1x Loop electrode 0°, Ø 27 Fr., monopolar | # 32-4310 |
| 1x Protection- and autoclaving tube for electrodes and cold knives | # 32-4230 |
| 1x Bladder syringe, 150 ml | # 32-4395-150 |
| 1x Ellik evacuator | # 32-4106 |
| 1x Monopolar cable, for use with monopolar Urology instruments, 3 m | # 11-1260UR |



24 / 26 Fr. Monopolar Resectoscopy Set

For Adults


- | | |
|--|---------------|
| 1x FullHD Cystoscope, Ø 4 mm, 30°, 300 mm | # 32-4020HD |
| 1x Resectoscopy sheath incl. obturator, continuous-flow, rotatable, Ø 24 / 26 Fr. | # 32-4205RR |
| 1x Visual obturator with Ø 8 Fr. working channel for Ø 24 / 26 Fr. resectoscopy sheath | # 32-4228 |
| 1x Monopolar working element, closed handle | # 32-4220 |
| 1x Loop electrode 30°, Ø 24 Fr., monopolar | # 32-4240 |
| 1x Ball electrode Ø 3 mm, Ø 24 Fr., monopolar | # 32-4250 |
| 1x Ball electrode Ø 5 mm, Ø 24 Fr., monopolar | # 32-4260 |
| 1x Knife electrode , Ø 24 Fr., monopolar | # 32-4265 |
| 1x Protection- and autoclaving tube for electrodes and cold knives | # 32-4230 |
| 1x Bladder syringe, 150 ml | # 32-4395-150 |
| 1x Ellik evacuator | # 32-4106 |
| 1x Monopolar cable, for use with monopolar Urology instruments, 3 m | # 11-1260UR |



19 / 22 Fr. Monopolar Resectoscopy Set

For Adolescents


- | | |
|---|---------------|
| 1x FullHD Cystoscope, Ø 2.9 mm, 30°, 300 mm | # 32-4022HD |
| 1x Resectoscopy sheath incl. obturator, continuous-flow, rotatable, Ø 19 / 22 Fr. | # 32-4292RR |
| 1x Visual obturator for Ø 19 / 22 Fr. resectoscopy sheath | # 32-4342 |
| 1x Monopolar working element, closed handle | # 32-4275 |
| 1x Roller barrel electrode Ø 3 mm, Ø 19 Fr., monopolar | # 32-4274 |
| 1x Loop electrode 30°, Ø 19 Fr., monopolar | # 32-4277 |
| 1x Knife electrode, Ø 19 Fr., monopolar | # 32-4279 |
| 1x Protection- and autoclaving tube for electrodes and cold knives | # 32-4230 |
| 1x Bladder syringe, 150 ml | # 32-4395-150 |
| 1x Ellik evacuator | # 32-4106 |
| 1x Monopolar cable, for use with monopolar Urology instruments, 3 m | # 11-1260UR |



23 Fr. Universal Laser Cystoscopy Set

Lateral Laser Emission for PVP / HoLAP / ThuVAP

- | | |
|---|-------------|
| 1x FullHD Cystoscope, Ø 2.9 mm, 30°, 300 mm | # 32-4002HD |
| 1x Laser cysto-urethroscope sheath incl. obturator, continuous-flow, Ø 23 Fr. | # 32-4600 |
| 1x Visual obturator for 32-4600 | # 32-4601 |
| 1x Telescope bridge with Ø 8 Fr. laser channel | # 32-4602 |




Compatible with:

Greenlight® Fibers MoXy with Ø 2.1 mm / max. Ø 2.3 mm, Lumenis® Xpeeda (Duotome) with Ø 2.5 mm, Quantana SideFiber® with Ø 2.22 mm (± 5%), Biolitec SideFiber 2100® with Ø 2.4 mm

24 / 26 Fr. Laser Resectoscopy Set

Lateral Laser Emission for PVP / HoLAP / ThuVAP

- | | |
|---|-------------|
| 1x FullHD Cystoscope, Ø 4 mm, 30°, 300 mm | # 32-4020HD |
| 1x Resectoscopy sheath incl. obturator, continuous-flow, rotatable, Ø 24 / 26 Fr. | # 32-4205RR |
| 1x Laser working insert, channel for laser for probes up to 7 Fr. | # 32-4610 |



Compatible with:

Greenlight® Fibers MoXy with Ø 2.1 mm / max. Ø 2.3 mm, Quantana SideFiber® with Ø 2.22 mm (± 5%), Biolitec® XCAVATOR with Glass Tip

24 / 26 Fr. Laser Resectoscopy Set, demountable

Straight Laser Emission for HoLEP / ThuLEP

- | | |
|--|----------------|
| 1x FullHD Cystoscope, Ø 4 mm, 30°, 300 mm | # 32-4020HD |
| 1x Resectoscopy sheath incl. obturator, continuous-flow, rotatable, Ø 24 / 26 Fr. | # 32-4205RR |
| 1x Adapter for insertion of laser probes | # 32-4607 |
| 1x Rubber reducer | # 11-1249-100 |
| 1x Laser working element, closed handle, passive, for exchangeable laser guiding pipes | # 32-4609 |
| 1x Guide pipe for laser probes Ø 0.8 mm, with retractor | # 32-4609-100R |
| 1x Guide pipe for laser probes Ø 1.2 mm, with retractor | # 32-4609-300R |
| 1x Guide pipe for laser probes Ø 1.4 mm, with retractor | # 32-4609-400R |
| 1x Guide pipe for laser probes Ø 0.8 mm | # 32-4609-100 |
| 1x Guide pipe for laser probes Ø 1.0 mm | # 32-4609-200 |
| 1x Guide pipe for laser probes Ø 1.2 mm | # 32-4609-300 |
| 1x Guide pipe for laser probes Ø 1.4 mm | # 32-4609-400 |




Compatible with:

Lumenis® Slimline, EZ, GI with max. Ø 1.2 mm (3.6 Fr.), Lisa Laser® with RigiFib-SV, Quanta® Laser with Cyber™, Boston Scientific® with AccuMax

24 / 26 Fr. Laser Resectoscopy Set, Ø 0.8 mm

Straight Laser Emission for HoLEP / ThuLEP

- | | |
|---|---------------|
| 1x FullHD Cystoscope, Ø 4 mm, 30°, 300 mm | # 32-4020HD |
| 1x Resectoscopy sheath incl. obturator, continuous-flow, rotatable, Ø 24 / 26 Fr. | # 32-4205RR |
| 1x Laser working element, closed handle, passive, channel for probes up to Ø 0.8 mm | # 32-4603 |
| 1x Adapter for insertion of laser probes | # 32-4607 |
| 1x Rubber reducer | # 11-1249-100 |




Compatible with:

Lumenis® Slimline, EZ, GI with max. Ø 1.2 mm (3.6 Fr.), Lisa Laser® with RigiFib-SV, Quanta® Laser with Cyber™, Boston Scientific® with AccuMax

24 / 26 Fr. Laser Resectoscopy Set, Ø 1.2 mm

Straight Laser Emission for HoLEP / ThuLEP

- | | |
|---|---------------|
| 1x FullHD Cystoscope, Ø 4 mm, 30°, 300 mm | # 32-4020HD |
| 1x Resectoscopy sheath incl. obturator, continuous-flow, rotatable, Ø 24 / 26 Fr. | # 32-4205RR |
| 1x Laser working element, closed handle, passive, channel for probes up to Ø 1.2 mm | # 32-4604 |
| 1x Adapter for insertion of laser probes | # 32-4607 |
| 1x Rubber reducer | # 11-1249-100 |



Compatible with:

Lumenis® Slimline, EZ, GI with max. Ø 1.2 mm (3.6 Fr.), Lisa Laser® with RigiFib-SV, Quanta® Laser with Cyber™, Boston Scientific® with AccuMax

19 / 22 Fr. Laser Resectoscopy Set, Ø 0.8 mm

Straight Laser Emission for HoLEP / ThuLEP

- | | |
|---|---------------|
| 1x FullHD Cystoscope, Ø 2.9 mm, 30°, 300 mm | # 32-4022HD |
| 1x Resectoscopy sheath incl. obturator, continuous-flow, rotatable, Ø 19 / 22 Fr. | # 32-4292RR |
| 1x Laser working element, closed handle, passive, channel for probes up to Ø 0.8 mm | # 32-4605 |
| 1x Adapter for insertion of laser probes | # 32-4607 |
| 1x Rubber reducer | # 11-1249-100 |




Compatible with:

Lumenis® Slimline, EZ, GI with max. Ø 1.2 mm (3.6 Fr.), Lisa Laser® with RigiFib-SV, Quanta® Laser with Cyber™, Boston Scientific® with AccuMax

19 / 22 Fr. Laser Resectoscopy Set, Ø 1.2 mm

Straight Laser Emission for HoLEP / ThuLEP

- | | |
|---|---------------|
| 1x FullHD Cystoscope, Ø 2.9 mm, 30°, 300 mm | # 32-4022HD |
| 1x Resectoscopy sheath incl. obturator, continuous-flow, rotatable, Ø 19 / 22 Fr. | # 32-4292RR |
| 1x Laser working element, closed handle, passive, channel for probes up to Ø 1.2 mm | # 32-4606 |
| 1x Adapter for insertion of laser probes | # 32-4607 |
| 1x Rubber reducer | # 11-1249-100 |




Compatible with:

Lumenis® Slimline, EZ, GI with max. Ø 1.2 mm (3.6 Fr.), Lisa Laser® with RigiFib-SV, Quanta® Laser with Cyber™, Boston Scientific® with AccuMax

17.5 / 18.5 Fr. Laser Resectoscopy Set, Ø 0.8 mm

Straight Laser Emission for HoLEP / ThuLEP

- | | |
|---|---------------|
| 1x FullHD Cystoscope, Ø 2.9 mm, 30°, 300 mm | # 32-4022HD |
| 1x Resectoscopy sheath incl. obturator, continuous-flow, rotatable, Ø 17.5 / 18.5 Fr. | # 32-4282RR |
| 1x Laser working element, closed handle, passive, channel for probes up to Ø 0.8 mm | # 32-4605 |
| 1x Adapter for insertion of laser probes | # 32-4607 |
| 1x Rubber reducer | # 11-1249-100 |




Compatible with:

Lumenis® Slimline, EZ, GI with max. Ø 1.2 mm (3.6 Fr.), Lisa Laser® with RigiFib-SV, Quanta® Laser with Cyber™, Boston Scientific® with AccuMax

17.5 / 18.5 Fr. Laser Resectoscopy Set, Ø 1.2 mm

Straight Laser Emission for HoLEP / ThuLEP

- | | |
|---|---------------|
| 1x FullHD Cystoscope, Ø 2.9 mm, 30°, 300 mm | # 32-4022HD |
| 1x Resectoscopy sheath incl. obturator, continuous-flow, rotatable, Ø 17.5 / 18.5 Fr. | # 32-4282RR |
| 1x Laser working element, closed handle, passive, channel for probes up to Ø 1.2 mm | # 32-4606 |
| 1x Adapter for insertion of laser probes | # 32-4607 |
| 1x Rubber reducer | # 11-1249-100 |




Compatible with:

Lumenis® Slimline, EZ, GI with max. Ø 1.2 mm (3.6 Fr.), Lisa Laser® with RigiFib-SV, Quanta® Laser with Cyber™, Boston Scientific® with AccuMax

Mechanical Lithotripsy Set


- | | |
|---|-------------|
| 1x FullHD Cystoscope, Ø 4 mm, 0°, 300 mm | # 32-4000HD |
| 1x FullHD Cystoscope, Ø 4 mm, 30°, 300 mm | # 32-4020HD |
| 1x Mauermayer stone punch, for small stones and fragments, QL | # 32-4069 |
| 1x Lithotripsy sheath incl. obturator, central valve, Ø 24 Fr., QL | # 32-4061 |
| 1x Visual obturator for lithotripsy sheath, Ø 24 Fr., Ø 7 Fr. working channel, QL | # 32-4062 |
| 1x Optical stone forceps, for small stones and stone fragments only | # 32-4066 |
| 1x Cysto-urethroscope sheath incl. obturator, 2 lateral stopcocks, Ø 25 Fr. | # 32-4090 |
| 1x Ellik evacuator | # 32-4106 |



21 Fr. Visual Urethrotomy Set

For Adults - NO HF!


- | | |
|--|---------------|
| 1x FullHD Cystoscope, Ø 4 mm, 0°, 300 mm | # 32-4000HD |
| 1x Urethrotomy sheath incl. obturator, 1 instrument channel, 2 lateral stopcocks, Ø 21 Fr. | # 32-4370 |
| 1x Monopolar working element, closed handle, passive | # 32-4220 |
| 1x Wave-cut cold knife, Ø 21 Fr. | # 32-4384 |
| 1x Half-round cold knife, Ø 21 Fr. | # 32-4389 |
| 1x Straight "Sachse" cold knife, Ø 21 Fr. | # 32-4390 |
| 1x Telescope bridge with 1 instrument channel, Ø 21 Fr. | # 32-4380 |
| 1x Urethrotome outer sheath for balloon catheter, use with 32-4370 | # 32-4375 |
| 1x Urethrotome outer sheath for continuous-flow, use with 32-4370 | # 32-4376 |
| 1x Bladder syringe, 150 ml | # 32-4395-150 |
| 1x Otis urethrotome | # 32-4391 |
| 1x Protection- and autoclaving tube for electrodes and cold knives | # 32-4230 |



15 Fr. Visual Urethrotomy Set

For Adolescents - NO HF!


- | | |
|--|--------------|
| 1x FullHD Cystoscope, Ø 2.9 mm, 0°, 300 mm | # 32-4002HD |
| 1x Urethrotomy sheath incl. obturator, 1 instrument channel, 2 lateral stopcocks, Ø 15 Fr. | # 32-4370-15 |
| 1x Working element, closed handle, passive, w/o HF-connector | # 32-4220UT |
| 1x Wave-cut cold knife, Ø 15 Fr. | # 32-4384-15 |
| 1x Half-round cold knife, Ø 15 Fr. | # 32-4389-15 |
| 1x Straight "Sachse" cold knife, Ø 15 Fr. | # 32-4390-15 |
| 1x Telescope bridge with 1 instrument channel with Ø 5 Fr., Ø 15 Fr. | # 32-4380-15 |



24 / 26 Fr. Short Nephroscopy Set

200 mm


- | | |
|---|--------------|
| 1x Short Nephroscope, Ø 4.8 mm, 6° viewing angle, 90° ocular, 200 mm | # 70-8201 |
| 1x Nephroscopy sheath incl. obturator, 1 stopcock, rotatable, Ø 24 Fr., 200 mm | # 70-8205-24 |
| 1x Nephroscopy sheath incl. obturator, 1 stopcock, rotatable, Ø 26 Fr., 200 mm | # 70-8205 |
| 1x Telescope dilation sheath set, 6 pcs., 9/12/15/18/21/24 Fr., 1x rigid and 1x flexible guide wire | # 70-8206 |
| 1x Dilator in 27 Fr. | # 70-8208 |
| 1x Dilator in 30 Fr. | # 70-8208-30 |
| 1x Stone grasping forceps, Ø 3 mm, rigid, ring handle, flushable, 400 mm | # 70-8220 |
| 1x Biopsy spoon, Ø 3 mm, rigid, ring handle, flushable, 400 mm | # 70-8225 |
| 1x Serrated grasping forceps, Ø 3 mm, rigid, ring handle, flushable, 400 mm | # 70-8230 |



24 / 26 Fr. Standard Nephroscopy Set

250 mm


- | | |
|---|---------------|
| 1x Standard Nephroscope, Ø 4.8 mm, 6° viewing angle, 90° ocular, 250 mm | # 70-8201-250 |
| 1x Nephroscopy sheath incl. obturator, 1 stopcock, rotatable, Ø 24 Fr., 250 mm | # 70-8204-24 |
| 1x Nephroscopy sheath incl. obturator, 1 stopcock, rotatable, Ø 26 Fr., 250 mm | # 70-8204 |
| 1x Telescope dilation sheath set, 6 pcs., 9/12/15/18/21/24 Fr., 1x rigid and 1x flexible guide wire | # 70-8206 |
| 1x Dilator in 27 Fr. | # 70-8208 |
| 1x Dilator in 30 Fr. | # 70-8208-30 |
| 1x Stone grasping forceps, Ø 3 mm, rigid, ring handle, flushable, 400 mm | # 70-8220 |
| 1x Biopsy spoon, Ø 3 mm, rigid, ring handle, flushable, 400 mm | # 70-8225 |
| 1x Serrated grasping forceps, Ø 3 mm, rigid, ring handle, flushable, 400 mm | # 70-8230 |



Compact Nephroscopy Set











250 mm

- | | |
|---|-------------------------|
| 1x Standard Nephroscope, Ø 4.8 mm, 6° viewing angle, 90° ocular, 200 mm | # 70-8200C |
| 1x Nephroscopy OP sheath incl. obturator, for use with Amplatz sheath, Ø 24 Fr. | # 70-8209
70-8209A |
| 1x Amplatz Sheath, for use with 70-8209 | |
| 1x Small stone grasping forceps, Ø 5 Fr., semi-rigid, ring handle, rotatable, 400 mm | # 70-8255SR |
| 1x Biopsy spoon, Ø 5 Fr., semi-rigid, ring handle, rotatable, 400 mm | # 70-8260SR |
| 1x Medium stone grasping forceps, Ø 5 Fr., semi-rigid, ring handle, rotatable, 400 mm | # 70-8275SR |




Mini Nephroscopy Set

220 mm

1x Mini Nephroscope, 220 mm	# 70-8200MINI	
1x Nephroscopy OP sheath incl. obturator, continuous-flow, Ø 15 / 18 Fr.	# 70-8211	
1x Nephroscopy OP sheath incl. obturator, continuous-flow, Ø 16.5 / 19.5 Fr.	# 70-8212	
1x Nephroscopy OP sheath incl. obturator, continuous-flow, Ø 21 / 24 Fr.	# 70-8213	
1x Dilator for nephroscopy OP sheath, 1 x cannulated, Ø 15 / 18 Fr.	# 70-8211D	
1x Dilator for nephroscopy OP sheath, 2 x cannulated, Ø 16.5 / 19.5 Fr.	# 70-8212D	
1x Dilator for nephroscopy OP sheath, 2 x cannulated, Ø 21 / 24 Fr.	# 70-8213D	
1x Small stone grasping forceps, Ø 5 Fr., semi-rigid, ring handle, rotatable, 400 mm	# 70-8255SR	
1x Biopsy spoon, Ø 5 Fr., semi-rigid, ring handle, rotatable, 400 mm	# 70-8260SR	
1x Medium stone grasping forceps, Ø 5 Fr., semi-rigid, ring handle, rotatable, 400 mm	# 70-8275SR	

FullHD Standard Uretero-Renoscapy Set

- | | |
|--|-----------|
| 1 x Ultra Thin HD Uretero-Renoscope, for instruments up to Ø 5 Fr., 310 mm | # 70-8001 |
| 1 x Dormia stone basket, Ø 5 Fr., 4 wires, 600 mm | # 70-8035 |
| 1 x Handle for dormia stone basket, Ø 5 Fr. | # 32-4211 |
| 1 x Small stone grasping forceps, Ø 5 Fr., flexible, 600 mm | # 70-8010 |
| 1 x Biopsy spoon, Ø 5 Fr., flexible, 600 mm | # 70-8015 |
| 1 x Medium stone grasping forceps, Ø 5 Fr., flexible, 600 mm | # 70-8020 |
| 1 x Large stone grasping forceps, Ø 5 Fr., flexible, 600 mm | # 70-8025 |



Thin HD Uretero-Renoscapy Set

- | | |
|---|-------------|
| 1x Ultra Thin HD Uretero-Renoscope, for instruments up to Ø 4 Fr., 310 mm | # 70-8001T |
| 1x Small stone grasping forceps, Ø 4 Fr., flexible, 600 mm | # 70-8010-4 |
| 1x Biopsy spoon, Ø 4 Fr., flexible, 600 mm | # 70-8015-4 |
| 1x Medium stone grasping forceps, Ø 4 Fr., flexible, 600 mm | # 70-8020-4 |
| 1x Large stone grasping forceps, Ø 4 Fr., flexible, 600 mm | # 70-8025-4 |



Ultra-Thin HD Uretero-Renoscropy Set

1x Ultra Thin HD Uretero-Renoscope, for instruments up to
Ø 3 Fr., 310 mm

70-8001UT



1x Small stone grasping forceps, Ø 3 Fr., flexible, 600 mm

70-8010-3



1x Medium stone grasping forceps, Ø 3 Fr., flexible, 600 mm

70-8020-3



1x Scissors, Ø 3 Fr., flexible, 600 mm

70-8030-3





Cystoscopy

Table of Contents

FullHD Cystoscopes / **p. 36 - 37**

Cysto-Urethroscope Sheaths & Obturators / **p. 38 - 40 & 52**

Optical Instruments, Telescope Bridges & Albarran Deflectors / **p. 41 - 43**

Flexible & Semi-Rigid Instruments / **p. 46 - 47 & 53 - 47**

Electrodes & Dormia Stone Baskets / **p. 48 - 49**

Compact Mini Cystoscopes / **p. 50 - 51**

Flexible Cystoscope / **p. 56 - 57**






FullHD Cystoscopes

Key Features

- Available in **FullHD resolution**
- **Sapphire glass lense** at proximal and distal end
- Nitrogen filled optical system to **avoid fogging** from the inside
- Three-tube construction for **more stability**
- Sealing by **laser welding** with **excellent** optical transmission
- **Autoclavable**



... more Features

- Made of high-quality stainless steel
- Sophisticated, anti-reflective optical system
- Powerful optical fibers
- Highly refractive lenses for optimum colour rendering even at high magnification
- Light cable adapters to Ackermann, Storz¹, Olympus¹, Wolf¹, ACMI¹
- High-quality glass cone for optimizing light input
- Standard C-Mount coupling
- Extremely precise centering of the optical and mechanical components for a brilliant and high-contrast image right up to the edge area the inside




FullHD Cystoscopes



	Diameter		Viewing angle				Item no.
	Ø 2.9 mm	Ø 4 mm	0°	12°	30°	70°	
FullHD Cystoscope, 300 mm	•		•				32-4002HD
FullHD Cystoscope, 300 mm		•	•				32-4000HD
FullHD Cystoscope, 300 mm	•			•			32-4015HD
FullHD Cystoscope, 300 mm		•		•			32-4010HD
FullHD Cystoscope, 300 mm	•				•		32-4022HD
FullHD Cystoscope, 300 mm		•			•		32-4020HD
FullHD Cystoscope, 300 mm		•				•	32-4030HD

Cysto-Urethroscope Sheaths incl. Obturator

2 Lateral Stopcocks / Adults & Adolescents



	Diameter	For use with scope		Color coding	Item no.
		Ø 2.9 mm	Ø 4 mm		
<p>Cysto-Urethroscope sheath incl. obturator, 2 lateral stopcocks</p> <p>includes Cysto-Urethroscope sheath, 2 lateral stopcocks, Ø 15 Fr. 32-4073S and Standard obturator for 32-4073S 32-40740</p>	Ø 15 Fr.	•		●	32-4073
<p>Cysto-Urethroscope sheath incl. obturator, 2 lateral stopcocks</p> <p>includes Cysto-Urethroscope sheath, 2 lateral stopcocks, Ø 17 Fr. 32-4070S and Standard obturator for 32-4070S 32-40700</p>	Ø 17 Fr.		•	●	32-4070
<p>Cysto-Urethroscope sheath incl. obturator, 2 lateral stopcocks</p> <p>includes Cysto-Urethroscope sheath, 2 lateral stopcocks, Ø 19 Fr. 32-4075S and Standard obturator for 32-4075S 32-40750</p>	Ø 19 Fr.		•	●	32-4075
<p>Cysto-Urethroscope sheath incl. obturator, 2 lateral stopcocks</p> <p>includes Cysto-Urethroscope sheath, 2 lateral stopcocks, Ø 21 Fr. 32-4080S and Standard obturator for 32-4080S 32-40800</p>	Ø 21 Fr.		•	●	32-4080
<p>Cysto-Urethroscope sheath incl. obturator, 2 lateral stopcocks</p> <p>includes Cysto-Urethroscope sheath, 2 lateral stopcocks, Ø 23 Fr. 32-4085S and Standard obturator for 32-4085S 32-40850</p>	Ø 23 Fr.		•	●	32-4085
<p>Cysto-Urethroscope sheath incl. obturator, 2 lateral stopcocks</p> <p>includes Cysto-Urethroscope sheath, 2 lateral stopcocks, Ø 25 Fr. 32-4090S and Standard obturator for 32-4090S 32-40900</p>	25 Fr.		•	○	32-4090

Cysto-Urethroscope Sheaths incl. Obturator

Central Valve / Adults & Adolescents



			Diameter	Color coding	Item no.
Cysto-Urethroscope sheath incl. obturator, central valve for use with Ø 4 mm scope			Ø 17 Fr.	●	32-4072
includes	Cysto-Urethroscope sheath, central valve, Ø 17 Fr.	32-4072S			
and	Standard obturator for 32-4072S	32-40700			
Cysto-Urethroscope sheath incl. obturator, central valve for use with Ø 4 mm scope			Ø 19 Fr.	●	32-4077
includes	Cysto-Urethroscope sheath, central valve, Ø 19 Fr.	32-4077S			
and	Standard obturator for 32-4077S	32-40750			
Cysto-Urethroscope sheath incl. obturator, central valve for use with Ø 4 mm scope			Ø 21 Fr.	●	32-4082
includes	Cysto-Urethroscope sheath, central valve, Ø 21 Fr.	32-4082S			
and	Standard obturator for 32-4082S	32-40800			
Cysto-Urethroscope sheath incl. obturator, central valve for use with Ø 4 mm scope			Ø 23 Fr.	●	32-4087
includes	Cysto-Urethroscope sheath, central valve, Ø 23 Fr.	32-4087S			
and	Standard obturator for 32-4087S	32-40850			
Cysto-Urethroscope sheath incl. obturator, central valve for use with Ø 4 mm scope			Ø 23 Fr.	●	32-4092
includes	Cysto-Urethroscope sheath, central valve, Ø 25 Fr.	32-4092S			
and	Standard obturator for 32-4092S	32-40900			

Standard Obturators for Cysto-Urethroscope Sheaths


For Use with Sheaths with 2 Lateral Stopcocks and Central Valve / Adults & Adolescents



	For use with scope		Diameter	Color coding	Item no.
	Ø 2.9 mm	Ø 4 mm			
Standard obturator for cysto-urethroscope sheath	•		15 Fr.	●	32-40740
Standard obturator for cysto-urethroscope sheath		•	17 Fr.	●	32-40700
Standard obturator for cysto-urethroscope sheath		•	19 Fr.	●	32-40750
Standard obturator for cysto-urethroscope sheath		•	21 Fr.	●	32-40800
Standard obturator for cysto-urethroscope sheath		•	23 Fr.	●	32-40850
Standard obturator for cysto-urethroscope sheath		•	25 Fr.	○	32-40900

Visual Obturators for Cysto-Urethroscope Sheaths

For Use with Sheaths with 2 Lateral Stopcocks and Central Valve / Adults & Adolescents



	Diameter	Color coding	Item no.
Visual obturator for cysto-urethroscope-sheath, for use with Ø 4 mm scope	17 Fr.	●	32-4071
Visual obturator for cysto-urethroscope-sheath, for use with Ø 4 mm scope	19 Fr.	●	32-4076
Visual obturator for cysto-urethroscope-sheath, for use with Ø 4 mm scope	21 Fr.	●	32-4081
Visual obturator for cysto-urethroscope-sheath, for use with Ø 4 mm scope	23 Fr.	●	32-4086
Visual obturator for cysto-urethroscope-sheath, for use with Ø 4 mm scope	25 Fr.	○	32-4091

Optical Grasping Forceps

For Use with 17 Fr. Cystoscopy Sheath or Resectoscopy Sheath plus 32-4067 / 32-4068



Item no.

Optical grasping forceps, 4.3 mm, double-action, for 30° scope

32-4050

Optical Biopsy Spoon

For Use with 17 Fr. Cystoscopy Sheath or Resectoscopy Sheath plus 32-4067 / 32-4068




Item no.

Optical biopsy spoon, 4.3 mm, double-action, spoon-type

32-4055

Optical Biopsy Spoon

For Use with 17 Fr. Cystoscopy Sheath or Resectoscopy Sheath plus 32-4067 / 32-4068



Item no.

Optical biopsy spoon, 2.3 mm, double-action, for 70° scope

32-4056

Optical Scissors

For Use with 17 Fr. Cystoscopy Sheath or Resectoscopy Sheath plus 32-4067 / 32-4068




Item no.

Optical scissors, hook-shaped, double-action

32-4060

Optical Biopsy Punch

For Use with 17 Fr. Cystoscopy Sheath or Resectoscopy Sheath plus 32-4067 / 32-4068



Item no.

Optical biopsy punch, through-cutting, for 30° scope

32-4065

Adapter for Optical Forceps



Item no.

Adapter for optical forceps 32-4066, use with Ackermann resectoscopy outer sheaths

32-4067

Adapter for Optical Forceps




Adapter for optical forceps, use with Ackermann resectoscopy inner sheaths

Item no.

32-4068

Telescope Bridge

For Diagnostics




Telescope bridge for diagnostics, suitable for 19 - 25 Fr. cysto-sheath

Item no.

32-4095

Telescope Bridge

With 1 Instrument Channel




Telescope bridge with 1 instrument channel, suitable for 19 - 25 Fr. cysto-sheath

Item no.

32-4100

Telescope Bridge

With 2 Instrument Channels



Item no.

Telescope bridge with 2 instrument channels, suitable for 19 - 25 Fr. cysto-sheath

32-4105

Albarran Deflector

With 1 Instrument Channel



Item no.

Albarran deflector with 1 instrument channel, lockable, suitable for cysto-sheaths from 17 Fr.

32-4110

Albarran Deflector

With 2 Instrument Channels



Item no.

Albarran deflector with 2 instrument channels, lockable, suitable for cysto-sheaths from 17 Fr.

32-4115


Instrument Compatibility Matrix

For Cysto-Urethroscope Sheaths with Instruments

Telescope bridge 32-4100, 1 instrument channel	Telescope bridge 32-4105, 2 instrument channels	Albarran 32-4110, 1 instrument channel	Albarran 32-4115, 2 instrument channels	Color coding of sheath
1 x 5 Fr.	2 x 4 Fr.	1 x 5 Fr.	2 x 4 Fr.	●
1 x 6 Fr.	2 x 5 Fr.	1 x 5 Fr.	2 x 4 Fr.	●
1 x 7 Fr.	2 x 5 Fr.	1 x 6 Fr.	2 x 5 Fr.	●
1 x 9 Fr.	2 x 6 Fr.	1 x 9 Fr.	2 x 5 Fr.	●
1 x max. 10 Fr.	2 x 8 Fr.	1 x max. 10 Fr.	2 x 8 Fr.	●

Biospy Forceps

For Adults & Adolescents



	Diameter			Type		Item no.
	Ø 5 Fr.	Ø 7 Fr.	Ø 9 Fr.	flexible	semi-rigid	
Biospy forceps, 400 mm	•			•		32-4121
Biospy forceps, 400 mm		•		•		32-4120
Biospy forceps, 400 mm			•	•		32-4120-9
Biospy forceps, 400 mm	•				•	32-4121SR
Biospy forceps, 400 mm		•			•	32-4120SR
Biospy forceps, 400 mm			•		•	32-4120-9SR

Foreign Body Forceps

For Adults & Adolescents



	Diameter			Type		Item no.
	Ø 5 Fr.	Ø 7 Fr.	Ø 9 Fr.	flexible	semi-rigid	
Foreign body forceps, 400 mm	•			•		32-4126
Foreign body forceps, 400 mm		•		•		32-4125
Foreign body forceps, 400 mm	•				•	32-4126SR
Foreign body forceps, 400 mm		•			•	32-4125SR

Serrated Grasping Forceps

For Adults & Adolescents



	Diameter			Type		Item no.
	Ø 5 Fr.	Ø 7 Fr.	Ø 9 Fr.	flexible	semi-rigid	
Serrated grasping forceps, 400 mm	•			•		32-4131
Serrated grasping forceps, 400 mm		•		•		32-4130
Serrated grasping forceps, 400 mm			•	•		32-4130-9
Serrated grasping forceps, 400 mm	•				•	32-4131SR
Serrated grasping forceps, 400 mm		•			•	32-4130SR
Serrated grasping forceps, 400 mm			•		•	32-4130-9SR

Scissors

For Adults & Adolescents



	Diameter			Type		Item no.
	Ø 5Fr.	Ø 7Fr.	Ø 9Fr.	flexible	semi-rigid	
Scissors, 400 mm	•			•		32-4136
Scissors, 400 mm		•		•		32-4135
Scissors, 400 mm			•	•		32-4135-9
Scissors, 400 mm	•				•	32-4136SR
Scissors, 400 mm		•		•	•	32-4135SR
Scissors, 400 mm			•		•	32-4135-9SR

Hook Scissors

For Adults & Adolescents



	Diameter		Type		Item no.
	Ø 5Fr.	Ø 9Fr.	flexible	semi-rigid	
Hook scissors, 400 mm	•		•		32-4140
Hook scissors, 400 mm		•	•		32-4140-9
Hook scissors, 400 mm	•			•	32-4140SR
Hook scissors, 400 mm		•		•	32-4140-9SR

DJ Grasping Forceps

For Adults & Adolescents



	Diameter			Type		Item no.
	Ø 5Fr.	Ø 7Fr.	Ø 9Fr.	flexible	semi-rigid	
DJ grasping forceps, 400 mm	•			•		32-4170
DJ grasping forceps, 400 mm		•		•		32-4170-7
DJ grasping forceps, 400 mm			•	•		32-4170-9
DJ grasping forceps, 400 mm	•				•	32-4170SR
DJ grasping forceps, 400 mm		•		•	•	32-4170-7SR
DJ grasping forceps, 400 mm			•		•	32-4170-9SR

Flexible Ball Electrodes, Monopolar

For Adults & Adolescents



	Diameter			Item no.
	Ø 3 Fr.	Ø 5 Fr.	Ø 7 Fr.	
Flexible ball electrode, monopolar, 470 mm	•			32-4202-3
Flexible ball electrode, monopolar, 470 mm		•		32-4202
Flexible ball electrode, monopolar, 470 mm			•	32-4150

Flexible Needle Electrodes, Monopolar

For Adults & Adolescents



	Diameter			Item no.
	Ø 3 Fr.	Ø 5 Fr.	Ø 7 Fr.	
Flexible needle electrode, monopolar, 470 mm	•			32-4201-3
Flexible needle electrode, monopolar, 470 mm		•		32-4201
Flexible needle electrode, monopolar, 470 mm			•	32-4155

Flexible Loop Electrodes, Monopolar

For Adults & Adolescents



	Diameter		Item no.
	Ø 5 Fr.	Ø 7 Fr.	
Flexible loop electrode, monopolar, 470 mm	•		32-4203-5
Flexible loop electrode, monopolar, 470 mm		•	32-4202

Flexible Ball Electrode, Bipolar

For Adults & Adolescents



Item no.

Flexible ball electrode, bipolar, Ø 5 Fr., 360 mm

32-4202BP

Flexible Needle Electrode, Bipolar

For Adults & Adolescents



Item no.

Flexible needle electrode, bipolar, Ø 5 Fr., 360 mm

32-4201BP

Dormia Stone Baskets

For Use with 32-4211 / 32-4212



	Diameter		Wires		Item no.
	Ø 5 Fr.	Ø 7 Fr.	4 wires	6 wires	
Dormia stone basket, 600 mm	•		•		32-4204
Dormia stone basket, 600 mm	•		•		32-4206
Dormia stone basket, 600 mm		•		•	32-4209
Dormia stone basket, 600 mm		•		•	32-4210

Handles for Dormia Stone Baskets



	Diameter		Item no.
	Ø 5 Fr.	Ø 7 Fr.	
Handle for dormia stone basket	•		32-4211
Handle for dormia stone basket		•	32-4212

Compact Mini Cystoscopes


Key Features

- **Soft-dilation** sheath in 3 sizes
- Working channel for use with **semi-rigid, flexible instruments** and high-flow cannula
- Large variety of all accessories for all kinds of **pediatric urology procedures**




... more Features

- Silicone valve for safe and easy handling of all accessories
- Light-weight engineering
- "Easy to use" concept
- 2 lateral suction-irrigation stop cocks




Compact Mini Cystoscopes



	Sheath diameter			Work. channel			Field of view		Length		Item no.
	Ø 4.9-6.5 Fr.	Ø 6-7 Fr.	Ø 7.5-9 Fr.	Ø 4 Fr.	Ø 5.1 Fr.	Ø 6.8 Fr.	80°	90°	150 mm	220 mm	
Compact Mini Cystoscope, 0°			•			•		•		•	32-4011
Compact Mini Cystoscope, 0°		•				•	•		•		32-4012
Compact Mini Cystoscope, 0°	•			•			•		•		32-4013

Cysto-Urethroscope Sheaths incl. Obturator

2 Lateral Stopcocks / Pediatrics



	Diameter	For use with scope		Color coding	Item no.
		Ø 2 mm	Ø 2.7 mm		
Cysto-Urethroscope sheath with 3 Fr. work. channel incl. obturator includes Cysto-Urethroscope sheath with working chan., Ø 9.5 Fr. 32-4511S and Standard obturator for 32-4511S 32-45110	Ø 9.5 Fr.	•		●	32-4511
Cysto-Urethroscope sheath with 3 Fr. work. channel incl. obturator includes Cysto-Urethroscope sheath with working chan., Ø 11 Fr. 32-4510S and Standard obturator for 32-4070S 32-45100	Ø 11 Fr.	•		●	32-4510
Cysto-Urethroscope sheath with 3 Fr. work. channel incl. obturator includes Cysto-Urethroscope sheath with working chan., Ø 13 Fr. 32-4298S and Standard obturator for 32-4075S 32-42980	Ø 13 Fr.		•	●	32-4298

Standard Obturators for Cysto-Urethroscope Sheaths


For Use with Sheaths with Working Channel / Pediatrics



	For use with scope		Diameter	Color coding	Item no.
	Ø 2 mm	Ø 2.7 mm			
Standard obturator for cysto-urethroscope-sheath	•		9.5 Fr.	●	32-45110
Standard obturator for cysto-urethroscope-sheath	•		11 Fr.	●	32-45100
Standard obturator for cysto-urethroscope-sheath		•	13 Fr.	●	32-42980

Biospy Forceps

For Pediatrics



	For use with scope			Diameter			Type		Item no.
	32-4013	32-4012	32-4011	Ø 3 Fr.	Ø 4 Fr.	Ø 5 Fr.	flexible	semi-rigid	
Biospy forceps, 280 mm	•			•			•		32-4162
Biospy forceps, 280 mm		•			•		•		32-4162-4
Biospy forceps, 280 mm			•			•	•		32-4162-5
Biospy forceps, 280 mm	•			•				•	32-4162SR
Biospy forceps, 280 mm		•			•		•		32-4162-4SR
Biospy forceps, 280 mm			•			•		•	32-4162-5SR

Foreign Body Forceps

for Pediatrics



	For use with scope		Diameter		Type		Item no.
	32-4012	32-4011	Ø 4 Fr.	Ø 5 Fr.	flexible	semi-rigid	
Foreign body forceps, 280 mm	•		•		•		32-4164-4
Foreign body forceps, 280 mm		•		•	•		32-4164-5
Foreign body forceps, 280 mm	•		•			•	32-4164-4SR
Foreign body forceps, 280 mm		•		•		•	32-4164-5SR

Serrated Grasping Forceps

For Pediatrics



	For use with scope			Diameter			Type		Item no.
	32-4013	32-4012	32-4011	Ø 3 Fr.	Ø 4 Fr.	Ø 5 Fr.	flexible	semi-rigid	
Serrated grasping forceps, 280 mm	•			•			•		32-4160
Serrated grasping forceps, 280 mm		•			•		•		32-4160-4
Serrated grasping forceps, 280 mm			•			•	•		32-4160-5
Serrated grasping forceps, 280 mm	•			•				•	32-4160SR
Serrated grasping forceps, 280 mm		•			•		•		32-4160-4SR
Serrated grasping forceps, 280 mm			•			•		•	32-4160-5SR

Scissors


For Pediatrics



	For use with scope		Diameter		Type		Item no.
	32-4012	32-4011	Ø 4 Fr.	Ø 5 Fr.	flexible	semi-rigid	
Scissors, 280 mm	•		•		•		32-4168-4
Scissors, 280 mm		•		•	•		32-4168-5
Scissors, 280 mm	•		•			•	32-4168-4SR
Scissors, 280 mm		•		•		•	32-4168-5SR

Flexible Ball Electrode


For Pediatrics



	Item no.
Flexible ball electrode, monopolar, Ø 3 Fr., 470 mm	32-4202-3

Flexible Needle Electrode

For Pediatrics



	Item no.
Flexible needle electrode, monopolar, Ø 3 Fr., 470 mm	32-4201-3



Flexible Cystoscope


Key Features

- **18.000 pixels** ensure highest image quality and
- A **100° viewing angle** allows for enhanced visibility
- The riveted deflection mechanism **reduces the risk of cracks** at the distal end
- Active dual angular deflection of **210° up and 120° down** for maximized versatility
- A **Ø 5.5 mm insertion tube** featuring a **Ø 2.3 mm working channel** allows for precise and delicate maneuverability




... more Features

- Total working length of 380 mm
- Light-weight design combined with a highly ergonomic handle for increased ease of use
- Compatible with all standard flexible instruments and camera systems
- Comes with purge valve, leakage tester and ETO cap



Flexible Cystoscope

For Use with Flexible Instruments



	Item no.
Flexible Cystoscope , Ø 16.5 Fr. sheath, Ø 6.9 Fr. work. chan., 210°/120° angulation, 100° field of view, 380 mm, includes suction valve, leakage tester and case	32-0000FLEX

Resectoscopy

Table of Contents

FullHD Cystoscopes / **p. 36 - 37**

Bipolar Resectoscopy / **p. 60 - 66**

Universal Resectoscopy / **p. 68 - 75**

Monopolar Resectoscopy / **p. 76 - 86**






Bipolar Resectoscopy

Key Features

- **True bipolar system** featuring a fully controlled flow of current eliminating the risk of unwanted thermal damage or electric exposure
- Specially insulated electrode tips ensure no **accidental tissue coagulation** by a second electrode
- **Convenient electrode replacement** without the need of removing the cable in situ
- **Open system** with enhanced compatibility properties




... more Features

- No risk of TUR syndrome by use of NaCl
- Single bipolar connector and quick assembly
- Excellent cutting, ignition and coagulation performance due to minimized distance between electrodes
- Premium, durable reusable electrodes



Bipolar Working Element, Passive


For Use with Ø 24 / 26 Fr. Bipolar Electrodes, Sheath and Ø 4 mm Scope



	Item no.
Bipolar working element , closed handle, passive	32-4220BP

Bipolar Working Element, Active


For Use with Ø 24 / 26 Fr. Bipolar Electrodes, Sheath and Ø 4 mm Scope



	Item no.
Bipolar working element , closed handle, active	32-4216BP

Bipolar Loop Electrode, 90°


For Use with 0° Scope



	Color coding	Item no.
Loop electrode, 90° , bipolar, Ø 24 / 26 Fr.	●	32-4235BP

Bipolar Loop Electrode, Retrograde


For Use with 30° Scope



	Color coding	Item no.
Loop electrode, retrograde, bipolar, Ø 24 / 26 Fr.	●	32-4240BP

Bipolar Loop Electrode, Straight


For Use with 0° / 12° / 30° / 70° Scope



	Color coding	Item no.
Loop electrode, straight, bipolar, Ø 24 / 26 Fr.	●	32-4245BP

Bipolar Ball Electrode

For Use with 0° / 12° / 30° / 70° Scope



	Color coding	Item no.
Ball electrode Ø 3 mm, bipolar, Ø 24 / 26 Fr.	●	32-4250BP

Bipolar Knife Electrode


For Use with 0° / 12° / 30° / 70° Scope



	Color coding	Item no.
Knife electrode, bipolar, Ø 24 / 26 Fr.	●	32-4265BP

Bipolar Roller Barrel Electrode

For Use with 0° / 12° / 30° / 70° Scope



	Color coding	Item no.
Roller barrel electrode, bipolar, Ø 24 / 26 Fr.	●	32-4266BP

Resectoscopy Sheath incl. Obturator

Continuous-flow, Rotatable / Adults



	Color coding	Item no.
Resectoscopy sheath incl. obturator, continuous-flow, rotatable, Ø 24 / 26 Fr., for use with Ø 4 mm scope	●	32-4205RR
includes Resectoscopy outer sheath, continuous-flow, rotatable		32-4205RR0S
and Resectoscopy inner sheath, continuous-flow, rotatable		32-4205RRIS
and Standard obturator for resectoscopy sheath		32-4227

Resectoscopy Sheath incl. Obturator

Continuous-flow / Adults



	Color coding	Item no.
Resectoscopy sheath incl. obturator , continuous-flow, Ø 24 / 26 Fr., for use with Ø 4 mm scope	●	32-4205
includes Resectoscopy outer sheath, continuous-flow		32-42050S
and Resectoscopy inner sheath, continuous-flow		32-4205IS
and Standard obturator for resectoscopy sheath		32-4227

Resectoscopy Sheath incl. Obturator

1 Stopcock / Adults



	Color coding	Item no.
Resectoscopy sheath incl. obturator , 1 stopcock, Ø 24 / 26 Fr., for use with Ø 4 mm scope	●	32-4222
includes Resectoscopy sheath, 1 stopcock		32-4222S
and Standard obturator for resectoscopy sheath		32-4227

Resectoscopy Sheath incl. Obturator

1 Stopcock / Adults



	Color coding	Item no.
Resectoscopy sheath incl. obturator , 1 stopcock, Ø 24 / 26 Fr., for use with Ø 4 mm scope	●	32-4205CV
includes Resectoscopy sheath, 1 stopcock		32-4205CVS
and Standard obturator for resectoscopy sheath		32-4227

Standard Resectoscopy Obturator



	Color coding	Item no.
Standard obturator for resectoscopy sheath, Ø 24 / 26 Fr.	●	32-4227

Visual Resectoscopy Obturator



	Color coding	Item no.
Visual obturator for resectoscopy sheath, Ø 24 / 26 Fr.	●	32-4228

Visual Resectoscopy Obturator, Working Channel



	Color coding	Item no.
Visual obturator for resectoscopy sheath with Ø 8 Fr. working channel, Ø 24 / 26 Fr.	●	32-4228-6

Deflecting Resectoscopy Obturator




	Color coding	Item no.
Deflecting obturator for resectoscopy sheath, Ø 24 / 26 Fr.	●	32-4226



Universal Resectoscopy


Key Features

- Combines **monopolar** and **bipolar Resectoscopy** in one system for complete versatility
- **Quick-Lock connector** for convenient assembly and disassembly
- Ergonomic, **light-weight handle design** with finger tang for maximum comfort
- **Complete range** of working elements and electrodes available




... more Features

- Monopolar mode by using a monopolar cable, neutral electrode and Glycine/ Sorbitol irrigation fluid
- Bipolar mode by using a bipolar cable and NaCl irrigation fluid



Universal Working Element, Passive


For Use with Ø 24 / 26 Fr. Universal Electrodes, Sheath and Ø 4 mm Scope



	Item no.
Universal working element , closed handle, passive	32-4220U

Universal Working Element, Active


For Use with Ø 24 / 26 Fr. Universal Electrodes, Sheath and Ø 4 mm Scope



	Item no.
Universal working element , closed handle, active	32-4216U

Universal Working Element, Passive

For Use with Ø 17.5 / 18.5 Fr. Universal Electrodes, Sheath and Ø 2.9 mm Scope



	Item no.
QL Universal working element , closed handle, passive	32-4275U

Universal Loop Electrode 90°


For Use with 0° Scope



	Color coding	Item no.
Loop electrode, straight, universal, Ø 24 / 26 Fr.		32-4235U

Universal Loop Electrodes, Retrograde


For Use with 30° Scope



	For use with sheath Ø 17.5 / 18.5 Fr. Ø 24 / 26 Fr.	Color coding	Item no.
Loop electrode, retrograde, universal	•		32-4287U
Loop electrode, retrograde, universal	•		32-4240U

Universal Loop Electrodes, Straight

For Use with 0° / 12° / 30° / 70° Scope



	For use with sheath Ø 17.5 / 18.5 Fr. Ø 24 / 26 Fr.	Color coding	Item no.
Loop electrode, straight, universal	•		32-4286U
Loop electrode, straight, universal	•		32-4245U

Universal Ball Electrodes


For Use with 0° / 12° / 30° / 70° Scope



	For use with sheath Ø 17.5 / 18.5 Fr. Ø 24 / 26 Fr.	Color coding	Item no.
Ball electrode , Ø 3 mm, universal	•	●	32-4288U
Ball electrode , Ø 3 mm, universal	•	●	32-4250U

Universal Knife Electrodes

For Use with 0° / 12° / 30° / 70° Scope



	For use with sheath Ø 17.5 / 18.5 Fr. Ø 24 / 26 Fr.	Color coding	Item no.
Knife electrode , universal	•	●	32-4289U
Knife electrode , universal	•	●	32-4265U

Universal Vaporization Electrode


For Use with 0° / 12° / 30° / 70° Scope



	Color coding	Item no.
Vaporization electrode , universal, Ø 24 / 26 Fr.	●	32-4255U

Universal Roller Barrel Electrodes

For Use with 0° / 12° / 30° / 70° Scope



	For use with sheath	Color coding	Item no.
	Ø 17.5 / 18.5 Fr. Ø 24 / 26 Fr.		
Roller Barrel electrode, Ø 3 mm, universal	•	●	32-4284U
Roller Barrel electrode, Ø 3 mm, universal	•	●	32-4266U

Resectoscopy Sheaths incl. Obturator

Continuous-flow, Rotatable / Adults & Adolescents



	Diameter	For use with scope	Color coding	Item no.
		Ø 2.9 mm Ø 4 mm		
Resectoscopy sheath incl. obturator, continuous-flow, rotatable includes Resectoscopy outer sheath, continuous-flow, rotatable 32-4205RR0S and Resectoscopy inner sheath, continuous-flow, rotatable 32-4205RRIS and Standard obturator for resectoscopy sheath 32-4227	Ø 24 / 26 Fr.	•	●	32-4205RR
QL Resectoscopy sheath incl. obturator, continuous-flow, rotatable includes Resectoscopy outer sheath, continuous-flow, rotatable 32-4282RR0S and Resectoscopy inner sheath, continuous-flow, rotatable 32-4282RRIS and Standard obturator for resectoscopy sheath 32-4344	Ø 17.5/18.5 Fr.	•	●	32-4282RR

Resectoscopy Sheaths incl. Obturator

Continuous-flow / Adults & Adolescents



			Diameter	For use with scope	Color coding	Item no.
				Ø 2.9 mm Ø 4 mm		
Resectoscopy sheath incl. obturator, continuous-flow			Ø 24 / 26 Fr.	•	●	32-4205
includes	Resectoscopy outer sheath, continuous-flow	32-42050S				
and	Resectoscopy inner sheath, continuous-flow	32-42051S				
and	Standard obturator for resectoscopy sheath	32-4227				
QL	Resectoscopy sheath incl. obturator, continuous-flow		Ø 17.5 / 18.5 Fr.	•	●	32-4281
includes	Resectoscopy outer sheath, continuous-flow	32-42810S				
and	Resectoscopy inner sheath, continuous-flow	32-42811S				
and	Standard obturator for resectoscopy sheath	32-4344				

Resectoscopy Sheaths incl. Obturator

1 Stopcock / Adults



			Color coding	Item no.
Resectoscopy sheath incl. obturator, 1 stopcock, Ø 24 / 26 Fr., for use with Ø 4 mm scope			●	32-4222
includes	Resectoscopy sheath, 1 stopcock	32-4222S		
and	Standard obturator for resectoscopy sheath	32-4227		

Resectoscopy Sheath incl. Obturator

1 Stopcock / Adults



		Color coding	Item no.
Resectoscopy sheath incl. obturator , 1 stopcock, Ø 24 / 26 Fr., for use with Ø 4 mm scope		●	32-4205CV
includes	Resectoscopy sheath, 1 stopcock		32-4205CVS
and	Standard obturator for resectoscopy sheath		32-4227

Standard Resectoscopy Obturators



	For use with sheath	Color coding	Item no.
Standard obturator for resectoscopy sheath	Ø 24 / 26 Fr.	●	32-4227
QL Standard obturator for resectoscopy sheath	Ø 17.5 / 18.5 Fr.	●	32-4344

Visual Resectoscopy Obturator




	Color coding	Item no.
Visual obturator for resectoscopy sheath, Ø 24 / 26 Fr.	●	32-4228

Visual Resectoscopy Obturator, Working Channel



	Diameter	Color coding	Item no.
Visual obturator for resectoscopy sheath with Ø 8 Fr. working channel	Ø 24 / 26 Fr.	●	32-4228-6
QL Visual obturator for resectoscopy sheath with Ø 5 Fr. working channel, QL	Ø 17.5 / 18.5 Fr.	●	32-4342-17

Deflecting Resectoscopy Obturator



	Color coding	Item no.
Deflecting obturator for resectoscopy sheath	●	32-4226

Monopolar Resectoscopy

Key Features

- Designed for the treatment of **adults, adolescents and children**
- Ergonomic, **light-weight handle design** with finger tang for maximum comfort
- **Open system** with enhanced compatibility properties
- **Complete range** of working elements and electrodes available




... more Features

- To be used with a monopolar cable, neutral electrode and Glycine/ Sorbitol irrigation fluid



Monopolar Working Element, Passive


For Use with Monopolar Electrodes, Sheath and Ø 2.9 mm / 4 mm Scope



	For use with electrodes			Color coding	Item no.
	Ø 19 Fr.	Ø 24 Fr.	Ø 27 Fr.		
Monopolar working element, closed handle, passive		•	•	● ●	32-4220
Monopolar working element, closed handle, passive	•			●	32-4275

Monopolar Working Element, Active

For Use with Monopolar Electrodes, Sheath and Ø 2.9 mm / 4 mm Scope



	For use with electrodes			Color coding	Item no.
	Ø 19 Fr.	Ø 24 Fr.	Ø 27 Fr.		
Monopolar working element, closed handle, active		•	•	● ●	32-4216
Monopolar working element, closed handle, active	•			●	32-4285

Monopolar Loop Electrodes, 90°

For Use with 0° Scope



	For use with sheath			Color coding	Item no.
	Ø 19 / 22 Fr.	Ø 24 / 26 Fr.	Ø 27 / 28.5 Fr.		
Loop electrode, 90°, monopolar			•	●	32-4295
Loop electrode, 90°, monopolar		•		●	32-4235
Loop electrode, 90°, monopolar	•			●	32-4235

Monopolar Loop Electrodes, Retrograde

For Use with 30° Scope



	For use with sheath			Color coding	Item no.
	Ø 19 / 22 Fr.	Ø 24 / 26 Fr.	Ø 27 / 28.5 Fr.		
Loop electrode, retrograde, monopolar			•	●	32-4300
Loop electrode, retrograde, monopolar		•		●	32-4240
Loop electrode, retrograde, monopolar	•			●	32-4277

Monopolar Loop Electrode, Straight

For Use with 0° / 12° / 30° / 70° Scope



	For use with sheath			Color coding	Item no.
	Ø 19 / 22 Fr.	Ø 24 / 26 Fr.	Ø 27 / 28.5 Fr.		
Loop electrode, straight, monopolar			•	●	32-4310
Loop electrode, straight, monopolar		•		●	32-4245
Loop electrode, straight, monopolar	•			●	32-4276

Monopolar Ball Electrodes


For Use with 0° / 12° / 30° / 70° Scope



	Ball diameter		For use with sheath		Color coding	Item no.
	Ø 3 mm	Ø 5 mm	Ø 19 / 22 Fr.	Ø 24 / 26 Fr.		
Ball electrode, monopolar	•			•	●	32-4250
Ball electrode, monopolar		•		•	●	32-4260
Ball electrode, monopolar	•		•		●	32-4278

Monopolar Knife Electrodes

For Use with 0° / 12° / 30° / 70° Scope



	For use with sheath	Color coding	Item no.
	Ø 19 / 22 Fr. Ø 24 / 26 Fr.		
Knife electrode, monopolar	•	•	32-4279
Knife electrode, monopolar	•	•	32-4265

Monopolar Vaporization Electrodes


For Use with 0° / 12° / 30° / 70° Scope



	Tip diameter	Color coding	Item no.
	Ø 3 mm Ø 5 mm		
Vaporization electrode, monopolar, for use with Ø 24 / 26 Fr. sheath	•	•	32-4255
Vaporization electrode, monopolar, for use with Ø 24 / 26 Fr. sheath	•	•	32-4256

Monopolar Cone Electrode

For Use with 0° / 12° / 30° / 70° Scope



	Color coding	Item no.
Cone electrode, monopolar, for use with Ø 24 / 26 Fr. sheath	•	32-4258

Monopolar Roller Barrel Electrodes

For Use with 0° / 12° / 30° / 70° Scope



	Ball diameter		For use with sheath		Color coding	Item no.
	Ø 3 mm	Ø 5 mm	Ø 19 / 22 Fr.	Ø 24 / 26 Fr.		
Roller Barrel electrode, monopolar	•			•	●	32-4266
Roller Barrel electrode, monopolar		•		•	●	32-4271
Roller Barrel electrode, monopolar	•		•		●	32-4274

Resectoscopy Sheaths incl. Obturator

Continuous-Flow, Rotatable / Adults & Adolescents



	Diameter	For use with scope		Color coding	Item no.
		Ø 2.9 mm	Ø 4 mm		
Resectoscopy sheath incl. obturator, continuous-flow, rotatable includes Resectoscopy outer sheath, continuous-flow, rotatable 32-4221RROS and Resectoscopy inner sheath, continuous-flow, rotatable 32-4221RRIS and Standard obturator for resectoscopy sheath 32-4231	Ø 27.5/28.5 Fr.		•	●	32-4221RR
Resectoscopy sheath incl. obturator, continuous-flow, rotatable includes Resectoscopy outer sheath, continuous-flow, rotatable 32-4205RROS and Resectoscopy inner sheath, continuous-flow, rotatable 32-4205RRIS and Standard obturator for resectoscopy sheath 32-4227	Ø 24 / 26 Fr.		•	●	32-4205RR
Resectoscopy sheath incl. obturator, continuous-flow, rotatable includes Resectoscopy outer sheath, continuous-flow, rotatable 32-4292RROS and Resectoscopy inner sheath, continuous-flow, rotatable 32-4292RRIS and Standard obturator for resectoscopy sheath 32-42910	Ø 19 / 22 Fr.	•		●	32-4292RR

Resectoscopy Sheaths incl. Obturator

Continuous-Flow / Adults & Adolescents



			Diameter	For use with scope	Color coding	Item no.
				Ø 2.9 mm Ø 4 mm		
Resectoscopy sheath incl. obturator, continuous-flow			Ø 27.5/28.5 Fr.	•	●	32-4221
includes	Resectoscopy outer sheath, continuous-flow	32-42210S				
and	Resectoscopy inner sheath, continuous-flow	32-42211S				
and	Standard obturator for resectoscopy sheath	32-4231				
Resectoscopy sheath incl. obturator, continuous-flow			Ø 24 / 26 Fr.	•	●	32-4205
includes	Resectoscopy outer sheath, continuous-flow	32-42050S				
and	Resectoscopy inner sheath, continuous-flow	32-42051S				
and	Standard obturator for resectoscopy sheath	32-4227				
Resectoscopy sheath incl. obturator, continuous-flow			Ø 19 / 22 Fr.	•	●	32-4290
includes	Resectoscopy outer sheath, continuous-flow	32-42900S				
and	Resectoscopy inner sheath, continuous-flow	32-42901S				
and	Standard obturator for resectoscopy sheath	32-42910				

Resectoscopy Sheaths incl. Obturator

1 Stopcock / Adults



			Diameter	Color coding	Item no.
Resectoscopy sheath incl. obturator, 1 stopcock, for use with Ø 4 mm scope			Ø 27.5/28.5 Fr.	●	32-4225
includes	Resectoscopy sheath, 1 stopcock	32-4225S			
and	Standard obturator for resectoscopy sheath	32-4231			
Resectoscopy sheath incl. obturator, 1 stopcock, for use with Ø 4 mm scope			Ø 24 / 26 Fr.	●	32-4222
includes	Resectoscopy sheath, 1 stopcock	32-4222S			
and	Standard obturator for resectoscopy sheath	32-4227			

Resectoscopy Sheaths incl. Obturator

1 Stopcock / Adults



			Diameter	Color coding	Item no.
Resectoscopy sheath incl. obturator, 1 stopcock			Ø 27.5 / 28.5 Fr.	●	32-4221CV
includes	Resectoscopy sheath, 1 stopcock	32-4221CVS			
and	Standard obturator for resectoscopy sheath	32-4231			
Resectoscopy sheath incl. obturator, 1 stopcock			Ø 24 / 26 Fr.	●	32-4205CV
includes	Resectoscopy sheath, 1 stopcock	32-4205CVS			
and	Standard obturator for resectoscopy sheath	32-4227			

Standard Resectoscopy Obturators



	For use with sheath	Color coding	Item no.
Standard obturator for resectoscopy sheath	Ø 27 / 28.5 Fr.	●	32-4231
Standard obturator for resectoscopy sheath	Ø 24 / 26 Fr.	●	32-4227
Standard obturator for resectoscopy sheath	Ø 19 / 22 Fr.	●	32-42910

Visual Resectoscopy Obturators




	For use with sheath	Color coding	Item no.
Visual obturator for resectoscopy sheath	Ø 27 / 28.5 Fr.	●	32-4240V0
Visual obturator for resectoscopy sheath	Ø 24 / 26 Fr.	●	32-4228
Visual obturator for resectoscopy sheath	Ø 19 / 22 Fr.	●	32-4342

Visual Resectoscopy Obturators, Working Channel



	Diameter	Color coding	Item no.
Visual obturator for resectoscopy sheath with Ø 9 Fr. working channel	Ø 27.5/28.5 Fr.	●	32-4240-9
Visual obturator for resectoscopy sheath with Ø 8 Fr. working channel	Ø 24 / 26 Fr.	●	32-4228-6


Deflecting Resectoscopy Obturators



	Diameter	Color coding	Item no.
Deflecting obturator for resectoscopy sheath	Ø 27.5/28.5 Fr.	●	32-4229
Deflecting obturator for resectoscopy sheath	Ø 24 / 26 Fr.	●	32-4226

Monopolar Working Element, Passive

For Use with Monopolar Electrodes, Sheath and Ø 2 mm /2.7 mm Scope



	For use with electrodes		Color coding	Item no.
	Ø 11 Fr.	Ø 13 Fr.		
Monopolar working element , closed handle, passive		•	●	32-4530
Monopolar working element , closed handle, passive	•		●	32-4296

Monopolar Loop Electrodes 90°, Large



	For use with sheath		Color coding	Item no.
	Ø 11 Fr.	Ø 13 Fr.		
Loop electrode 90° , large version, monopolar	•		●	32-4311
Loop electrode 90° , large version, monopolar		•	●	32-4301

Monopolar Loop Electrodes 90°, Small




	For use with sheath		Color coding	Item no.
	Ø 11 Fr.	Ø 13 Fr.		
Loop electrode 90° , small version, monopolar	•		●	32-4312
Loop electrode 90° , small version, monopolar		•	●	32-4302

Monopolar Ball Electrodes



	For use with sheath		Color coding	Item no.
	Ø 11 Fr.	Ø 13 Fr.		
Ball electrode 90°, monopolar	•		●	32-4313
Ball electrode 90°, monopolar		•	●	32-4303

Monopolar Knife Electrodes



	For use with sheath		Color coding	Item no.
	Ø 11 Fr.	Ø 13 Fr.		
Knife electrode, monopolar	•		●	32-4314
Knife electrode, monopolar		•	●	32-4304

Monopolar Hook Electrodes



	For use with sheath		Color coding	Item no.
	Ø 11 Fr.	Ø 13 Fr.		
Hook electrode, monopolar	•		●	32-4313WOB
Hook electrode, monopolar		•	●	32-4303WOB

Resectoscopy Sheaths incl. Obturator

1 Stopcock / Pediatrics



			Diameter	For use with scope	Color coding	Item no.
				Ø 2 mm Ø 2.7 mm		
Resectoscopy sheath incl. obturator, 1 stopcock			Ø 11 Fr.	•	●	32-4520
includes	Resectoscopy sheath, 1 stopcock	32-4520S				
and	Standard obturator for resectoscopy sheath	32-45200				
Resectoscopy sheath incl. obturator, 1 stopcock			Ø 13 Fr.	•	●	32-4297
includes	Resectoscopy sheath, 1 stopcock	32-4297S				
and	Standard obturator for resectoscopy sheath	32-42970				

Resectoscopy Obturator

Pediatrics



	For use with sheath	Color coding	Item no.
	Ø 11 Fr. Ø 13 Fr.		
Standard obturator for resectoscopy sheath	•	●	32-45200
Standard obturator for resectoscopy sheath	•	●	32-42970



Lithotripsy

Table of Contents

FullHD Cystoscopes / **p. 36 - 37**

Lithotripsy Sheath & Obturators / **p. 90**

Lithotripsy Speciality Instruments / **p. 91 - 92**





Lithotripsy Sheath incl. Obturator

Central Valve, Quick Lock



	Diameter	For use with scope	Color coding	Item no.
		Ø 4 mm, 0° Ø 4 mm, 30°		
Lithotripsy sheath incl. obturator, central valve, Quick Lock	Ø 24 Fr.	• •	●	32-4061
includes	Lithotripsy sheath, central valve, Quick-Lock			32-4061S
and	Standard obturator, Quick Lock, for 32-4061S			32-4063

Lithotripsy Standard Obturator

For Use with 32-4061S



	Item no.
Standard obturator for lithotripsy sheath, Ø 24 Fr., Quick-Lock	32-4063

Lithotripsy Visual Obturator


For Use with 32-4061S



	Item no.
Visual obturator for lithotripsy sheath, Ø 24 Fr., , Ø 7 Fr. working channel, Quick-Lock	32-4062

Mauermayer Stone Punch

For Use with 32-4061 / 32-4062 / 32-4063




Item no.

Mauermayer stone punch, for small stones and fragments, Quick-Lock

32-4069

Stone Lithotripter "The Beast"



Item no.

Stone lithotripter, demountable into 3 parts, Quick-Lock

32-4064

Optical Stone Forceps

For Use with 25 Fr. Cystoscopy Sheath plus 32-4067



Item no.

Optical stone forceps, for small stones and stone fragments only

32-4066

Adapter for Optical Stone Forceps



Adapter for optical stoneforceps 32-4066, use with Ackermann resectoscopy outer sheaths

Item no.
32-4067



Urethrotomy

Table of Contents

FullHD Cystoscopes / **p. 36 - 37**

Working Elements / **p. 96 & 101**

Urethrotomy Sheaths & Obturators / **p. 97 - 98 & 101**

Cold Knives / **p. 99 - 100 & 102 - 103**


Otis Urethrotome / **p. 100**





Monopolar Working Element, Passive

For Use with Cold Knives, without HF !



	Item no.
Monopolar working element , closed handle, passive, for use with Ø 21 Fr. cold knives without HF	32-4220

Monopolar Working Element, Active

For Use with Cold Knives, without HF !



	Item no.
Monopolar working element , closed handle, active, for use with Ø 21 Fr. cold knives without HF	32-4216

Working Element, Passive

For Use with Cold Knives, without HF !



	Item no.
Working element without HF connection , closed handle, passive, for use with Ø 15 Fr. cold knives	32-4220UT

Urethrotomy Sheaths incl. Obturator

2 Lateral Stopcocks / Adults & Adolescents



			Diameter		For use with scope		Item no.
			Ø 15 Fr.	Ø 21 Fr.	Ø 2.9 mm	Ø 4 mm	
Urethrotomy sheath incl. obturator, 2 lateral stopcocks				•		•	32-4370
includes	Urethrotomy sheath, 2 lateral stopcocks	32-4370S					
and	Standard obturator for 32-4370S	32-43700					
Urethrotomy sheath incl. obturator, 2 lateral stopcocks			•		•		32-4370-15
includes	Urethrotomy sheath, 2 lateral stopcocks	32-4370-15S					
and	Standard obturator for 32-4370-15S	32-4370-150					

Urethrotomy Outer Sheath for Balloon Catheter



	Item no.
Urethrotome outer sheath for balloon catheter, for use with 32-4370	32-4375

Urethrotomy Outer Sheath for Continuous-Flow



	Item no.
Urethrotome outer sheath for balloon catheter, for use with 32-4370	32-4376


Standard Urethrotomy Obturator

for Adults / Adolescents



	Diameter		For use with scope		Item no.
	Ø 15 Fr.	Ø 21 Fr.	Ø 2.9 mm	Ø 4 mm	
Standard obturator for urethrotomy sheath		•		•	32-43700
Standard obturator for urethrotomy sheath	•		•		32-4370-150

Visual Urethrotomy Obturator



	Item no.
Visual obturator for urethrotomy sheath, for use with 32-4370	32-4372

Wave-cut Cold Knives



	Diameter		For use with scope		Item no.
	Ø 15 Fr.	Ø 21 Fr.	Ø 2.9 mm	Ø 4 mm	
Wave-cut cold knife for visual urethrotome		•		•	32-4384
Wave-cut cold knife for visual urethrotome	•		•		32-4384-15

Half-round Cold Knife



	Item no.
Half-round cold knife for visual urethrotome, Ø 21 Fr., for use with Ø 4 mm scope	32-4389

Straight "Sachse" Cold Knives



	Diameter		For use with scope		Item no.
	Ø 15 Fr.	Ø 21 Fr.	Ø 2.9 mm	Ø 4 mm	
Straight "Sachse" cold knife for visual urethrotome		•		•	32-4390
Straight "Sachse" cold knife for visual urethrotome	•		•		32-4390-15


Hook-shaped Cold Knives



	Diameter		For use with scope		Item no.
	Ø 15 Fr.	Ø 21 Fr.	Ø 2.9 mm	Ø 4 mm	
Hook-shaped cold knife for visual urethrotome		•		•	32-4388
Hook-shaped cold knife for visual urethrotome	•		•		32-4388-15

Telescope Bridge

With 1 Instrument Channel



	Diameter		For use with scope		Item no.
	Ø 15 Fr.	Ø 21 Fr.	Ø 2.9 mm	Ø 4 mm	
Telescope bridge with 1 instrument channel		•		•	32-4380
Telescope bridge with 1 instrument channel	•		•		32-4380-15


Otis Urethrotome



	Item no.
Otis urethrotome	32-4391

Monopolar Working Elements, Passive


For Use with Cold Knives, without HF !



	For use with electrodes		Color coding	Item no.
	Ø 11 Fr.	Ø 13 Fr.		
Monopolar working element , closed handle, passive		•	●	32-4530
Monopolar working element , closed handle, passive	•		●	32-4296

Urethrotomy Sheaths, 2 Lateral Stopcocks

For Use with 32-4370-130 / 32-4370-110 / Pediatrics



	Diameter		For use with scope		Color coding	Item no.
	Ø 11 Fr.	Ø 13 Fr.	Ø 2 mm	Ø 2.7 mm		
Urethrotomy sheath incl. obturator , 2 lateral stopcocks		•		•	●	32-4370-13
includes Urethrotomy sheath, 2 lateral stopcocks 32-4370-13S						
and Standard obturator for 32-4370-13S 32-4370-130						
Urethrotomy sheath incl. obturator , 2 lateral stopcocks	•		•		●	32-4370-11
includes Urethrotomy sheath, 2 lateral stopcocks 32-4370-11S						
and Standard obturator for 32-4370-11S 32-4370-110						

Urethrotomy Obturators

For Use with 32-4370-13S / 32-4370-11S / Pediatrics



	Diameter		For use with scope		Color coding	Item no.
	Ø 11 Fr.	Ø 13 Fr.	Ø 2 mm	Ø 2.7 mm		
Standard obturator for urethrotomy sheath		•		•	• (Red)	32-4370-130
Standard obturator for urethrotomy sheath	•		•		• (Green)	32-4370-110

Half-round Cold Knives



	Diameter		For use with scope		Item no.
	Ø 11 Fr.	Ø 13 Fr.	Ø 2 mm, 187 mm	Ø 2.7 mm, 187 mm	
Half-round cold knife for visual urethrotome		•		•	32-4389-13
Half-round cold knife for visual urethrotome	•		•		32-4389-11

Sickle-shaped Cold Knives



	Diameter		For use with scope		Item no.
	Ø 11 Fr.	Ø 13 Fr.	Ø 2 mm, 187 mm	Ø 2.7 mm, 187 mm	
Sickle-shaped cold knife for visual urethrotome		•		•	32-4387-13
Sickle-shaped cold knife for visual urethrotome	•		•		32-4387-11


Hook-shaped Cold Knives



	Diameter		For use with scope		Item no.
	Ø 11 Fr.	Ø 13 Fr.	Ø 2 mm, 187 mm	Ø 2.7 mm, 187 mm	
Hook-shaped cold knife for visual urethrotome		•		•	32-4388-13
Hook-shaped cold knife for visual urethrotome	•		•		32-4388-11

Telescope Bridge

For Diagnostics



	Diameter		For use with scope		Item no.
	Ø 11 Fr.	Ø 13 Fr.	Ø 2 mm, 187 mm	Ø 2.7 mm, 187 mm	
Telescope bridge for pediatric visual urethrotome		•		•	32-4380-13
Telescope bridge for pediatric visual urethrotome	•		•		32-4380-11

Nephroscopy

Table of Contents

Short & Standard Nephroscopy / **p. 106 - 111**

Compact Nephroscopy / **p. 112 - 115**

Mini Nephroscopy / **p. 116 - 119**





Short & Standard Nephroscopes

Key Features

- **Slim and robust** design with \varnothing 4.8 mm working channel
- For use with **rigid Nephroscopy instruments** in 400 or 500 mm length
- Available with **45° or 90° angled ocular** and 200 mm or 250 mm working length




... more Features

- Suitable sheaths incl. obturator with a diameter of 24 / 26 Fr.
- Viewing angle of 6°




Short & Standard Nephroscopes, 45°



	Length		Item no.
	200 mm	250 mm	
Short Nephroscope , Ø 4.8 mm, 6° viewing angle, 45° ocular	•		70-8200
Standard Nephroscope , Ø 4.8 mm, 6° viewing angle, 45° ocular		•	70-8200-250

Short & Standard Nephroscopes, 90°



	Length		Item no.
	200 mm	250 mm	
Short Nephroscope , Ø 4.8 mm, 6° viewing angle, 90° ocular	•		70-8201
Standard Nephroscope , Ø 4.8 mm, 6° viewing angle, 90° ocular		•	70-8201-250

Nephroscopy Sheaths incl. Obturator

1 Lateral Stopcock, Rotatable / Short & Standard Nephroscopy



	Diameter		Use with scope		Color coding	Item no.
	Ø 24 Fr.	Ø 26 Fr.	in 200 mm	in 250 mm		
Nephroscopy sheath incl. obturator, 1 stopcock, rotatable includes Nephroscopy sheath, 1 stopcock, rotatable 70-8204S-24 and Standard obturator for 70-8204S-24 70-82040-24	•			•	●	70-8204-24
Nephroscopy sheath incl. obturator, 1 stopcock, rotatable includes Nephroscopy sheath, 1 stopcock, rotatable 70-8204S and Standard obturator for 70-8204S 70-82040		•		•	●	70-8204
Nephroscopy sheath incl. obturator, 1 stopcock, rotatable includes Nephroscopy sheath, 1 stopcock, rotatable 70-8205S-24 and Standard obturator for 70-8205S-24 70-82050-24	•			•	●	70-8205-24
Nephroscopy sheath incl. obturator, 1 stopcock, rotatable includes Nephroscopy sheath, 1 stopcock, rotatable 70-8205S and Standard obturator for 70-8205S 70-82050		•		•	●	70-8205

Nephroscopy Obturators


For Use with Sheath with 1 Lateral Stopcock, Rotatable / Short & Standard Nephroscopy



	Diameter		Color coding	Item no.
	Ø 24 Fr.	Ø 26 Fr.		
Nephroscopy obturator, cannulated	•		●	70-82040-24
Nephroscopy obturator, cannulated		•	●	70-82040
Nephroscopy obturator, cannulated	•		●	70-82050-24
Nephroscopy obturator, cannulated		•	●	70-82050

Stone Grasper

Short & Standard Nephroscopy



Stone grasping forceps, Ø 3 mm, rigid, ring handle, flushable, 400 mm

Item no.
70-8220

Biopsy Spoon

Short & Standard Nephroscopy



Biopsy spoon, Ø 3 mm, rigid, ring handle, flushable, 400 mm

Item no.
70-8225

Serrated Grasper

Short & Standard Nephroscopy



Serrated grasping forceps, Ø 3 mm, rigid, ring handle, flushable, 400 mm

Item no.
70-8230

Scissors

Short & Standard Nephroscopy



Scissors, Ø 3 mm, rigid, ring handle, flushable, 400 mm

Item no.

70-8235

Fenestrated Stone Grasper

Short & Standard Nephroscopy




Fenestrated stone grasping forceps, Ø 3 mm, spring handle, 400 mm

Item no.

70-8210

3-Prong Grasper

Short & Standard Nephroscopy



3-Prong grasping forceps, Ø 3 mm, spring handle, 400 mm

Item no.

70-8215

Straight Knife

Short & Standard Nephroscopy



Straight knife, Ø 3 mm, 3-ring handle, 400 mm	Item no. 70-8240
--	-----------------------------------

Sickle Knife

Short & Standard Nephroscopy



Sickle knife, Ø 3 mm, 3-ring handle, 400 mm	Item no. 70-8245
--	-----------------------------------

Suction Tube

Short & Standard Nephroscopy



Suction tube, Ø 3 mm, 400 mm	Item no. 70-8250
-------------------------------------	-----------------------------------

Compact Nephroscope


Key Features

- **Continuous-flow**-system
- Ø 18.8 Fr. sheath diameter and integrated rod lens system for **finest image quality**
- **Ideal for bariatric applications** due to 220 mm working length



... more Features

- Slim and robust design
- Oval sheath and round working channel
- Automatic valve
- Viewing angle of 6° allows field of view of 80°



Compact Nephroscope, 45°



	Item no.
Compact Nephroscope 45° , Ø 18.8 Fr., 12.9 Fr. work. ch., 6° viewing angle, 80° field of view, 220 mm	70-8200C

Nephroscopy OP Sheaths incl. Obturator

For Use with Scope 70-8200C & 70-8209A / **Compact Nephroscopy**



	Color coding	Item no.
Nephroscopy OP sheath incl. obturator, for use with Amplatz sheath, Ø 24 Fr.	●	70-8209

Amplatz Sheath for Nephroscopes


For Use with Scope 70-8200C & 70-8209A / **Compact Nephroscopy**



	Color coding	Item no.
Amplatz sheath , cannulated, Ø 27 Fr.	●	70-8209A

Stone Grasping Forceps

For Mini / Compact Nephroscopy




Stone grasping forceps for small stones, Ø 5 Fr., semi-rigid, ring handle, rotatable, 400 mm

Item no.

70-8255SR

Biopsy Spoon

For Mini / Compact Nephroscopy



Biopsy spoon, Ø 5 Fr., semi-rigid, ring handle, rotatable, 400 mm

Item no.

70-8260SR

Scissors

For Mini / Compact Nephroscopy



Scissors, Ø 5 Fr., semi-rigid, ring handle, rotatable, 400 mm

Item no.

70-8265SR

Stone Grasping Forceps

For Mini / Compact Nephroscopy




Stone grasping forceps for medium stones, Ø 5 Fr., semi-rigid, ring handle, rotatable, 400 mm

Item no.

70-8275SR

Stone Grasping Forceps

For Mini / Compact Nephroscopy



Stone grasping forceps for large stones, Ø 5 Fr., semi-rigid, ring handle, rotatable, 400 mm

Item no.

70-8280SR

Mini Nephroscope

Key Features

- Extremely slim design - 12 Fr. Mini Nephroscope **minimizes kidney trauma**
- **FullHD 40.000 pixel** resolution
- **Optic creates a flushing stream** which removes stone fragments
- **Continuous-flow**-system




... more Features

- Suitable operating sheaths have diameters of only Ø 15 Fr., Ø 16.5 Fr. or Ø 21 Fr.
- Working channel with Ø 6.6 Fr. for instruments up to Ø 5 Fr.
- Direction of view of 12° and field of view of 80°
- Kidney stones are fragmented with a lithotripter or laser lithotripsy system




Mini Nephroscope



	Item no.
Mini Nephroscope, Ø 12 Fr. sheath, Ø 6.6 Fr. work. chan., 12° viewing angle, 80° field of view, 220 mm	70-8200MINI

Nephroscopy OP Sheaths incl. Obturator

For Use with Scope 70-8200MINI & 70-8211D / 70-8212D / 70-8213D / Mini Nephroscopy



	Diameter			Color coding	Item no.
	Ø 15/18 Fr.	Ø 16.5/19.5 Fr.	Ø 21/24 Fr.		
Nephroscopy OP sheath incl. obturator, continuous-flow	•			•	70-8213
Nephroscopy OP sheath incl. obturator, continuous-flow		•		•	70-8212
Nephroscopy OP sheath incl. obturator, continuous-flow			•	•	70-8211

Dilators for Nephroscopy OP Sheaths

For Use with Sheaths above / Mini Nephroscopy



	Cannulated		Diameter			Color coding	Item no.
	1x	2x	Ø 15/18 Fr.	Ø 16.5/19.5 Fr.	Ø 21/24 Fr.		
Dilator for nephroscopy OP sheath, cannulated	•		•			•	70-8213D
Dilator for nephroscopy OP sheath, cannulated		•		•		•	70-8212D
Dilator for nephroscopy OP sheath, cannulated		•			•	•	70-8211D

Sealing Applicator for Nephroscopy OP Sheath


For Use with Sheaths above / **Mini Nephroscopy**



	Color coding	Item no.
Applicator for sealing , Ø 15/18 Fr., 2 parts in combination, sheath and plunger	●	70-8214

Stone Grasping Forceps


For Mini / Compact Nephroscopy



	Item no.
Stone grasping forceps for small stones , Ø 5 Fr., semi-rigid, ring handle, rotatable, 400 mm	70-8255SR

Biopsy Spoon

For Mini / Compact Nephroscopy



	Item no.
Biopsy spoon , Ø 5 Fr., semi-rigid, ring handle, rotatable, 400 mm	70-8260SR

Scissors

For Mini / Compact Nephroscopy



Scissors, Ø 5 Fr., semi-rigid, ring handle, rotatable, 400 mm

Item no.

70-8265SR

Stone Grasping Forceps

For Mini / Compact Nephroscopy



Stone grasping forceps for medium stones, Ø 5 Fr., semi-rigid, ring handle, rotatable, 400 mm

Item no.

70-8275SR

Stone Grasping Forceps

For Mini / Compact Nephroscopy



Stone grasping forceps for large stones, Ø 5 Fr., semi-rigid, ring handle, rotatable, 400 mm

Item no.

70-8280SR

Dilators

For Mini / Compact Nephroscopy



	Diameter							Item no.	
	Ø 9 Fr.	Ø 12 Fr.	Ø 15 Fr.	Ø 18 Fr.	Ø 21 Fr.	Ø 24 Fr.	Ø 27 Fr.		Ø 30 Fr.
Dilator	•							70-8206-09	
Dilator		•						70-8206-12	
Dilator			•					70-8206-15	
Dilator				•				70-8206-18	
Dilator					•			70-8206-21	
Dilator						•		70-8206-24	
Dilator							•	70-8208	
Dilator								•	70-8208-30

Guide Wires



	Type		Item no.
	flexible	rigid	
Guide wire for telescope dilation	•		70-8206-FF
Guide wire for telescope dilation		•	70-8206-FS

Telescope Dilation Sheath Set



			Item no.
Telescope dilation sheath set			70-8206
includes	Dilator in Ø 9 Fr.	70-8206-09	
and	Dilator in Ø 12 Fr.	70-8206-12	
and	Dilator in Ø 15 Fr.	70-8206-15	
and	Dilator in Ø 18 Fr.	70-8206-18	
and	Dilator in Ø 21 Fr.	70-8206-21	
and	Dilator in Ø 24 Fr.	70-8206-24	
and	Guide wire, flexible	70-8206-FF	
and	Guide wire, rigid	70-8206-FS	

Uretero-Renoscropy


Table of Contents

Semi-flexible Uretero-Renoscopes / **p. 124 - 125**

Semi-rigid and Flexible Instruments / **p. 126 - 129**

Flexible Uretero-Renoscope / **p. 130 - 131**






Uretero-Renoscopes

Key Features

- Available in **FullHD and HD resolution**
- **Extra-slim** working channel
- **Soft tapered sheath** allows for easy introduction and **atraumatic dilatation**
- **Semi-flexible** sheath




... more Features

- Available in different working lengths, suitable for upper and lower URS
- 45° angled light cable connector and instrument channels ensure ergonomic handling
- Improved dilatation and continuous suction enhance the field of vision during bleeding
- Leakage-free working channels due to a silicone valve with 2 parallel sealings



Standard FullHD Uretero-Renoscopes



	Sheath diameter	Length	Item no.
	Ø 7.5-9-12-13.2 Fr. Ø 7.5-9-12 Fr.	310 mm 425 mm	
FullHD Standard Uretero-Renoscope, 0°, Ø 6.8 Fr. work. ch., 90° field of view	•	•	70-8000
FullHD Standard Uretero-Renoscope, 0°, Ø 6.8 Fr. work. ch., 90° field of view	•	•	70-8001

Thin HD Uretero-Renoscopes



	Sheath diameter	Length	Item no.
	Ø 6-7.5-9-10.2 Fr. Ø 6-7-9 Fr.	310 mm 425 mm	
HD Thin Uretero-Renoscope, 0°, Ø 5.1 Fr. work. ch., 80° field of view	•	•	70-8000T
HD Thin Uretero-Renoscope, 0°, Ø 5.1 Fr. work. ch., 80° field of view	•	•	70-8001T

Ultra Thin HD Uretero-Renoscopes



	Sheath diameter	Length	Item no.
	Ø 4.9-6.5-8.1-9.3 Fr. Ø 4.9-6.5-8.1 Fr.	310 mm 425 mm	
HD Ultra Thin Uretero-Renoscope, 0°, Ø 4 Fr. work. ch., 80° field of view	•	•	70-8000UT
HD Ultra Thin Uretero-Renoscope, 0°, Ø 4 Fr. work. ch., 80° field of view	•	•	70-8001UT

Stone Grasping Forceps

For Small Stones



	Diameter			Type		Item no.
	Ø 3 Fr.	Ø 4 Fr.	Ø 5 Fr.	flexible	semi-rigid	
Stone grasping forceps, 600 mm	•			•		70-8010-3
Stone grasping forceps, 600 mm		•		•		70-8010-4
Stone grasping forceps, 600 mm			•	•		70-8010
Stone grasping forceps, 600 mm	•				•	70-8010-3SR
Stone grasping forceps, 600 mm		•		•	•	70-8010-4SR
Stone grasping forceps, 600 mm			•		•	70-8010SR

Stone Grasping Forceps


For Medium Stones



	Diameter		Type		Item no.
	Ø 4 Fr.	Ø 5 Fr.	flexible	semi-rigid	
Stone grasping forceps, 600 mm	•		•		70-8015-4
Stone grasping forceps, 600 mm		•	•		70-8015
Stone grasping forceps, 600 mm	•			•	70-8015-4SR
Stone grasping forceps, 600 mm		•		•	70-8015SR


Stone Grasping Forceps

For Large Stones




	Diameter		Type		Item no.
	Ø 4 Fr.	Ø 5 Fr.	flexible	semi-rigid	
Stone grasping forceps, 600 mm	•		•		70-8025-4
Stone grasping forceps, 600 mm		•	•		70-8025
Stone grasping forceps, 600 mm	•			•	70-8025-4SR
Stone grasping forceps, 600 mm		•		•	70-8025SR

Biopsy Spoon



	Diameter			Type		Item no.
	Ø 3 Fr.	Ø 4 Fr.	Ø 5 Fr.	flexible	semi-rigid	
Biopsy spoon, 600 mm	•			•		70-8020-3
Biopsy spoon, 600 mm		•		•		70-8020-4
Biopsy spoon, 600 mm			•	•		70-8020
Biopsy spoon, 600 mm	•				•	70-8020-3SR
Biopsy spoon, 600 mm		•		•	•	70-8020-4SR
Biopsy spoon, 600 mm			•		•	70-8020SR

Scissors



	Diameter			Type		Item no.
	Ø 3 Fr.	Ø 4 Fr.	Ø 5 Fr.	flexible	semi-rigid	
Scissors, 600 mm	•			•		70-8030-3
Scissors, 600 mm		•		•		70-8030-4
Scissors, 600 mm			•	•		70-8030
Scissors, 600 mm	•				•	70-8030-3SR
Scissors, 600 mm		•		•	•	70-8030-4SR
Scissors, 600 mm			•		•	70-8030SR

Dormia Stone Basket



	Item no.
Dormia stone basket, Ø 5 Fr., 4 wires, 600 mm	70-8035

Handle for Dormia Stone Basket



Handle for dormia stone basket, Ø 5 Fr.

Item no.

32-4211

Cleaning Brush



Cleaning brush, 600 mm

Item no.

70-8040

URS Adapter

For Use with URS Scopes



URS adapter

Item no.

70-8000-100

Silicone Cap Luer Lock

For Use with URS Scopes



	Item no. 70-8000A
Silicone cap luer lock	

Flexible Uretero-Renoscope


Key Features

- Enhanced active dual angular deflection of **275° up** and **275° down** for maximized versatility
- **Memory effect** in sheath for convenient readjustment
- Large **viewing angle of 88°** for increased visibility
- Total **working length of 680 mm** featuring a **Ø 1.2 mm working channel**



...more Features

- Distal tip diameter of Ø 2.5mm / 7.5 Fr.
- Ø 2.8 mm / 8.4 Fr. insertion tube diameter
- Light-weight design combined with a highly ergonomic handle for increased ease of use
- Compatible with all standard flexible instruments and camera systems
- Comes with leakage tester and ETO cap



Flexible Uretero-Renoscope

for use with flexible instruments



	Item no.
<p>Flexible Uretero-Renoscope, Ø 8.4 Fr. sheath, Ø 3.6 Fr. work. chan., 275° angulation, 60° field of view, 680 mm, includes suction valve, leakage tester and case</p>	70-0000FLEX

Urology Accessories

Table of Contents

Speciality Instruments / **p. 134 - 136**

Connection Cables / **p. 137 - 138**





Bougies by Dittel, Straight

More Sizes available in the ROLAN General Surgery Catalogue / Chapter 30



	Diameter	Item no.
Dittel Bougie , straight, without channel	Ø 9 Fr.	32-4400
Dittel Bougie , straight, without channel	Ø 10 Fr.	32-4401
Dittel Bougie , straight, without channel	Ø 12 Fr.	32-4402
Dittel Bougie , straight, without channel	Ø 13 Fr.	32-4404
Dittel Bougie , straight, without channel	Ø 15 Fr.	32-4405
Dittel Bougie , straight, without channel	Ø 16 Fr.	32-4406
Dittel Bougie , straight, without channel	Ø 18 Fr.	32-4407
Dittel Bougie , straight, without channel	Ø 19 Fr.	32-4408
Dittel Bougie , straight, without channel	Ø 21 Fr.	32-4409
Dittel Bougie , straight, without channel	Ø 22 Fr.	32-4410
Dittel Bougie , straight, without channel	Ø 24 Fr.	32-4411
Dittel Bougie , straight, without channel	Ø 25 Fr.	32-4412
Dittel Bougie , straight, without channel	Ø 27 Fr.	32-4413
Dittel Bougie , straight, without channel	Ø 28 Fr.	32-4414
Dittel Bougie , straight, without channel	Ø 30 Fr.	32-4415


Bougies by Dittel, Curved

More Sizes available in the ROLAN General Surgery Catalogue / Chapter 30




	Diameter	Item no.
Dittel Bougie , curved, without channel	ø 8 Fr.	32-4429
Dittel Bougie , curved, without channel	ø 9 Fr.	32-4430
Dittel Bougie , curved, without channel	ø 10 Fr.	32-4431
Dittel Bougie , curved, without channel	ø 12 Fr.	32-4432
Dittel Bougie , curved, without channel	ø 13 Fr.	32-4433
Dittel Bougie , curved, without channel	ø 15 Fr.	32-4434
Dittel Bougie , curved, without channel	ø 16 Fr.	32-4435
Dittel Bougie , curved, without channel	ø 18 Fr.	32-4436
Dittel Bougie , curved, without channel	ø 19 Fr.	32-4437
Dittel Bougie , curved, without channel	ø 21 Fr.	32-4438
Dittel Bougie , curved, without channel	ø 22 Fr.	32-4439
Dittel Bougie , curved, without channel	ø 24 Fr.	32-4440
Dittel Bougie , curved, without channel	ø 25.5 Fr.	32-4441
Dittel Bougie , curved, without channel	ø 27 Fr.	32-4442
Dittel Bougie , curved, without channel	ø 28.5 Fr.	32-4443
Dittel Bougie , curved, without channel	ø 30 Fr.	32-4444

Bladder Syringe



	Item no.
Bladder syringe, 150 ml	32-4395-150

Strauss Penis Clamp



	Item no.
Strauss penis clamp	32-4420





Ellik Evacuator



	Item no.
Ellik evacuator with rubber ball and locking	32-4106
includes Glass retainer only for Ellik evacuator	32-4106E
and Rubber ball for Ellik evacuator	32-4106-100
and Tubing with locking adapter for Ellik evacuator	32-4106-200





Bipolar Resectoscopy Cables

For Use with Bipolar Resectoscope

		Generator side			Item no.
		2-Banana	Storz ¹ / Erbe ¹	Martin ¹ / Berchtold ¹ / Aesculap ¹ GK55 and GK 60	
Round plug instrument side 	3 m length				
		•			13-1641UR
			•		13-1643UR
				•	13-1644UR






Universal Resectoscopy Cables

For Use with Universal Resectoscopes

		Generator side			Item no.
		2-Banana	Storz ¹ / Erbe ¹	Martin ¹ / Berchtold ¹ / Aesculap ¹ GK55 and GK 60	
Round plugs, 2 instrument sides 	3 m length				
		•			13-1641U
			•		
				•	13-1644U

Monopolar Resectoscopy Cables

For Use with Monopolar Resectoscopes

		Generator side				Item no.
		Ackermann / Erbe T ¹ / Martin ¹ / Berchthold ¹	Erbe ACC ¹ / ICC-VIO ¹	Bovie ¹ / Valleylab ¹ / Conmed ¹	Martin ¹ / Berchtold ¹ / Aesculap ¹	
						
		3 m length	3 m length	3 m length	3 m length	
Ø 2 mm instrument side, round plug 	•					11-1260UR
		•				11-1261UR
				•		11-1262UR
					•	11-1263UR



Alphanumeric Index

Item Numbers at a Glance

Numeric Item #			
11-1249-100	16, 17, 18, 19,	32-4068	43
	20, 21, 22	32-4069	23, 91, 93
11-1260UR	6, 11, 12, 13, 138	32-4070	6, 38
11-1261UR	138	32-40700	38, 39, 40
11-1262UR	138	32-4070S	38, 52
11-1263UR	138	32-4071	40
13-1641U	9, 10, 137	32-4072	39
13-1641UR	8, 137	32-4072S	38, 39
13-1643U	137	32-4073	38
13-1643UR	137	32-4073S	38
13-1644U	137	32-40740	38, 40
13-1644UR	137	32-4075	38
32-0000FLEX	57	32-40750	38, 39, 40
32-4000HD	23, 24, 37, 93	32-4075S	38, 52
32-4002HD	25, 37	32-4076	40
32-4010HD	37	32-4077	39
32-4011	51, 53, 54	32-4077S	39
32-4012	7, 51, 53, 54	32-4080	38
32-4013	51, 53	32-40800	38, 39, 40
32-4015HD	37	32-4080S	38
32-4020HD	6, 8, 9, 11, 12,	32-4081	40
	16, 17, 18, 23, 37, 93	32-4082	39
32-4022HD	10, 13, 19, 20,	32-4082S	39
	21, 22, 37	32-4085	6, 38
32-4030HD	37	32-40850	38, 39, 40
32-4050	41	32-4085S	38
32-4055	6, 41	32-4086	40
32-4056	41	32-4087	39
32-4060	42	32-4087S	39
32-4061	23, 90, 93	32-4090	23, 38, 93
32-4061S	90	32-40900	38, 39, 40
32-4062	23, 90, 93	32-4090S	38
32-4063	90	32-4091	40
32-4064	91	32-4092	39
32-4065	42	32-4092S	39
32-4066	23, 91, 93	32-4095	43
32-4067	42, 91, 92	32-4100	43, 45
		32-4105	6, 44, 45
		32-4106	6, 11, 12, 13, 23, 93,
			136
		32-4106-100	136
		32-4106-200	136
		32-4106E	136
		32-4110	44, 45
		32-4115	6, 44, 45
		32-4120	6, 46
		32-4120-9	46
		32-4120-9SR	46
		32-4120SR	46
		32-4121	46
		32-4121SR	46
		32-4125	46
		32-4125SR	46
		32-4126	46
		32-4126SR	46
		32-4130	6, 46
		32-4130-9	46
		32-4130-9SR	46
		32-4130SR	46
		32-4131	46
		32-4131SR	46
		32-4135	47
		32-4135-9	47
		32-4135-9SR	47
		32-4135SR	47
		32-4136	6, 47
		32-4136SR	47
		32-4140	47
		32-4140-9	47
		32-4140-9SR	47
		32-4140SR	47
		32-4150	6, 48
		32-4155	48
		32-4160	7, 53
		32-4160-4	53
		32-4160-4SR	53

32-4160-5	53	32-4212	49	32-4255U	71
32-4160-5SR	53	32-4216	77, 96	32-4256	79
32-4160SR	53	32-4216BP	61	32-4258	79
32-4162	7, 53	32-4216U	69	32-4260	12, 78
32-4162-4	53	32-4220	11, 12, 24, 77, 96	32-4265	12, 79
32-4162-4SR	53	32-4220BP	8, 61	32-4265BP	8, 63
32-4162-5	53	32-4220U	9, 69	32-4265U	9, 71
32-4162-5SR	53	32-4220UT	25, 96	32-4266	80
32-4162SR	53	32-4221	81	32-4266BP	63
32-4164-4	53	32-4221CV	82	32-4266U	72
32-4164-4SR	53	32-4221CVS	82	32-4271	80
32-4164-5	53	32-4221IS	81	32-4274	13, 80
32-4164-5SR	53	32-4221OS	81	32-4275	13, 77
32-4168-4	54	32-4221RR	11, 80	32-4275U	10, 69
32-4168-4SR	54	32-4221RRIS	80	32-4276	78
32-4168-5	54	32-4221RROS	80	32-4277	13, 78
32-4168-5SR	54	32-4222	64, 73, 81	32-4278	78
32-4170	47	32-4222S	64, 73, 81	32-4279	13, 79
32-4170-7	47	32-4225	81	32-4281	73
32-4170-7SR	47	32-4225S	81	32-4281OS	73
32-4170-9	47	32-4226	66, 75, 83	32-4281IS	73
32-4170-9SR	47	32-4227	63, 64, 65, 72, 73,	32-4282RR	10, 21, 22, 72
32-4170SR	47		74, 80, 81, 82	32-4282RRIS	72
32-4201	48	32-4228	12, 65, 74, 83	32-4282RROS	72
32-4201-3	48, 54	32-4228-6	8, 9, 65, 75, 83	32-4284U	72
32-4201BP	49	32-4229	83	32-4285	77
32-4202	6, 48	32-4230	8, 9, 10, 11, 12, 13,	32-4286U	70
32-4202-3	48, 54		24	32-4287U	10, 70
32-4202BP	49	32-4231	80, 81, 82	32-4288U	10, 71
32-4203-5	48	32-4235	77	32-4289U	10, 71
32-4204	6, 49	32-4235BP	61	32-4290	81
32-4205	64, 73, 81	32-4235U	70	32-4290IS	81
32-4205CV	64, 74, 82	32-4240	12, 78	32-4290OS	81
32-4205CVS	64, 74, 82	32-4240-9	11, 83	32-42910	80, 81, 82
32-4205IS	64, 73, 81	32-4240BP	8, 62	32-4292IS	81
32-4205OS	64, 73, 81	32-4240U	9, 70	32-4292RR	13, 19, 20, 80
32-4205RR	8, 9, 12, 16, 17,	32-4240VO	83	32-4292RRIS	80
	18, 63, 72, 80	32-4245	78	32-4292RROS	80
32-4205RRIS	63, 72, 80	32-4245BP	62	32-4295	11, 77
32-4205RROS	63, 72, 80	32-4245U	70	32-4296	84, 101
32-4206	49	32-4250	12, 78	32-4297	86
32-4209	49	32-4250BP	8, 62	32-42970	86
32-4210	49	32-4250U	9, 71	32-4297S	86
32-4211	6, 30, 49, 128	32-4255	79	32-4298	52

32-42980	52	32-4389-13	102	32-4520	86
32-4298S	52	32-4389-15	25	32-45200	86
32-4300	11, 78	32-4390	24, 99	32-4520S	86
32-4301	84	32-4390-15	25, 99	32-4530	84, 101
32-4302	84	32-4391	24, 100	32-4603	17
32-4303	85	32-4395-150	6, 11, 12, 13, 24, 136	32-4604	18
32-4303WOB	85			32-4605	19, 21
32-4304	85	32-4400	134	32-4606	20, 22
32-4310	11, 78	32-4401	134	32-4607	16, 17, 18, 19, 20, 21, 22
32-4312	84	32-4402	134		
32-4313	85	32-4404	134	32-4609	16
32-4313WOB	85	32-4405	134	32-4609-100	16
32-4314	85	32-4406	134	32-4609-100R	16
32-4342	13, 83	32-4407	134	32-4609-200	16
32-4342-17	10, 75	32-4408	134	32-4609-300	16
32-4344	72, 73, 74	32-4409	134	32-4609-300R	16
32-4370	24, 97, 98	32-4410	134	32-4609-400	16
32-4370-11	101	32-4411	134	32-4609-400R	16
32-4370-110	101, 102	32-4412	134	32-43011	84
32-4370-11S	101	32-4413	134	32-43014	85
32-4370-13	101	32-4414	134	70-0000FLEX	131
32-4370-130	101, 102	32-4415	134	70-8000	125
32-4370-13S	101	32-4420	136	70-8000-100	128
32-4370-15	25, 97	32-4429	135	70-8000A	129
32-4370-150	97, 98	32-4430	135	70-8000T	125
32-4370-15S	97	32-4431	135	70-8000UT	125
32-43700	97, 98	32-4432	135	70-8001	30, 124, 125
32-4370S	97	32-4433	135	70-8001T	31, 125
32-4372	98	32-4434	135	70-8001UT	32, 125
32-4375	24, 97	32-4435	135	70-8010	30, 126
32-4376	24, 98	32-4436	135	70-8010-3	32, 126
32-4380	24, 100	32-4437	135	70-8010-3SR	126
32-4380-11	103	32-4438	135	70-8010-4	31, 126
32-4380-13	103	32-4439	135	70-8010-4SR	126
32-4380-15	25, 100	32-4440	135	70-8010SR	126
32-4384	24, 99	32-4441	135	70-8015	30, 126
32-4384-15	25, 99	32-4442	135	70-8015-4	31, 126
32-4387-11	102	32-4443	135	70-8015-4SR	126
32-4387-13	102	32-4444	135	70-8015SR	126
32-4388	100	32-4510	52	70-8020	30, 127
32-4388-11	103	32-45100	52	70-8020-3	32, 127
32-4388-13	103	32-4510S	52	70-8020-3SR	127
32-4388-15	100	32-4511	52	70-8020-4	31, 127
32-4389	24, 99	32-45110	52	70-8020-4SR	127
32-4389-11	102	32-4511S	52	70-8020SR	127

70-8025	30, 126	70-8208	26, 27, 120
70-8025-4	31, 126	70-8208-30	26, 27, 120
70-8025-4SR	126	70-8209	28, 113
70-8025SR	126	70-8209A	28, 113
70-8030	127	70-8210	110
70-8030-3	32, 127	70-8211	29, 117
70-8030-3SR	127	70-8211D	29, 117
70-8030-4	127	70-8212	29, 117
70-8030-4SR	127	70-8212D	29, 117
70-8030SR	127	70-8213	29, 117
70-8035	30, 127	70-8213D	29, 117
70-8040	128	70-8214	118
70-8200	107	70-8215	110
70-8200-250	107	70-8220	26, 27, 109
70-8200C	28, 113	70-8225	26, 27, 109
70-8200MINI	29, 117	70-8230	26, 27, 109
70-8201	26, 107	70-8235	110
70-8201-250	27, 107	70-8240	111
70-8204	27, 108	70-8245	111
70-8204-24	27, 108	70-8250	111
70-82040	108	70-8255SR	28, 29, 114, 118
70-82040-24	108	70-8260SR	28, 29, 114, 118
70-8204S	108	70-8265SR	114, 119
70-8204S-24	108	70-8275SR	28, 29, 115, 119
70-8205	26, 108	70-8280SR	115, 119
70-8205-24	26, 108		
70-82050	108		
70-82050-24	108		
70-8205S	108		
70-8205S-24	108		
70-8206	26, 27, 121		
70-8206-09	120, 121		
70-8206-12	120, 121		
70-8206-15	120, 121		
70-8206-18	120, 121		
70-8206-21	120, 121		
70-8206-24	120, 121		
70-8206-FF	120, 121		
70-8206-FS	120, 121		

¹Any companies mentioned in this catalogue including their brand names and related products are protected under intellectual property laws. The use of them does not imply any affiliation with or endorsement by them.

²Specifications, design and accessories in this catalogue are subject to change without any notice or obligation on the part of the manufacturer.

Copyright by Ackermann Instrumente GmbH. All rights reserved. Specifications, design and accessories in this catalogue are subject to change without any notice or obligation on the part of the manufacturer. No part may be copied or reproduced in any form. Revision 10/21



www.ackermannsurgical.com



Ackermann Instrumente GmbH
Eisenbahnstrasse 65-67
78604 Rietheim-Weilheim, Germany