

electrosurgery

a compact solution for private practice



ES120

Electrosurgical device

ES120 is an electrosurgical device which uses high-frequency current flow for cutting and coagulating skin lesions.

Versatility of application

The unit allows to perform procedures in monopolar and bipolar modes. It excels in small surgical, gynecological, stomatological field and many others.

Automatic power regulation

The device features fully automatic adjustment of power output based on real-time monitoring of working parameters. Desired effect on tissue is repeatedly consistent in every working conditions.

Comfort

Small size and aesthetic design of ES120 along with its reliability and easy to use features ensure the highest comfort of work.

Safety

Start-up AutoTest of generator and plugged-in accessories as well as NEM system ensures maximum safety both for an operator and a patient.

Output configuration


ES120 is equipped with two 6-pin (SDS) outputs or standard 3-pin monopolar and 2-pin bipolar outputs.

Areas of applications:

- small surgical procedures
- oncology (skin lesions removal, samples for tests)
- gynecology
- laryngology
- dentistry










Advantages of ES120:

- small size and aesthetic design
- efficiency
- precision
- low operating costs



accessories for ES120


CE
1011

| | | |
|--|---------|---|
|  | 520-500 | Knife electrode straight, 25 x 3.5mm, handle 4mm |
|  | 521-420 | Ball electrode straight, 4mm, handle 4mm |
|  | 521-220 | Wire loop electrode, 10mm, handle 4mm |
|  | 524-020 | Ball electrode straight, 4mm, L: 100mm, handle 4mm |
|  | 525-000 | Loop electrode, 20 x 15mm, L: 120mm, handle 4mm |
|  | 526-400 | Triangle loop electrode, 10 x 25mm, L: 115mm, handle 4mm |
|  | 322-145 | Electrode handle 4mm, large, 2 switches, cable 4m, 3-pin plug |
|  | 605-014 | Bipolar forceps, angled, 1mm, L: 160 mm |
|  | 351-031 | Bipolar cable, L: 3m, straight connector, 2-pin 29mm plug |

non-stick




380-030 Disposable neutral electrode cable,
380-050 L: 3m or 5m, flat plug



812-80H Disposable neutral electrode EMED SAFE,
hydrogel, split, for adults and children
176x122mm, 110 cm², 10x5pcs.



100-301 One-pedal footswitch, for ES120,
cable 5m, 6-pin plug



100-302 2-button footswitch,
cable 5m, 6-pin plug

Full range of products are available in our catalogue of accessories.
See www.emed.pl for details.

REF NO. 2017.10.K19.EN

EMED SP. Z O. O. SP. K. | ul. Ryżowa 69a | 05-816 Opacz Kolonia | Poland |
tel: +48 (22) 723 08 00 | fax: +48 (22) 723 00 81 | email: export@emed.pl | www.emed.pl |